

Seeing **beyond**

ZEISS

ZEISS ClearMind Lenses
See clearly.
Think clearly.



This image was created with Generative AI.

Vision is our most dominant sense.

An estimated **80%** of the sensory information we rely on in everyday environments is **visual**. This makes clear vision important for living in our fast-paced, information-packed world.

Our eyes are bombarded with visual stimuli, like billboards, pop-ups, digital imagery, screens and notifications every day, many of these forcing us to multitask.

If your vision is unclear, your brain is forced to work harder. Blurred vision **can increase cognitive load** because the brain has to *fill in the gaps*.

This can reduce your concentration and the ability to stay focused, ultimately impacting your overall well-being.

At ZEISS, we realise that clear vision is more than just eyesight — it's a vital aspect of people's quality of life. This is why we dedicated our efforts to designing **our most sophisticated optical lenses for your eyes and mind — ZEISS ClearMind lenses.**

ZEISS ClearMind lenses

are the first of their kind from ZEISS. They provide extremely clear vision, helping to reduce your cognitive load and promote your overall well-being.*

* Opinion-based survey among wearers of ZEISS ClearMind lenses in Germany, Italy, India & China (N=298, Single Vision, Digital & Progressive lenses), 2025, Carl Zeiss Vision International GmbH, DE (unpublished, data on file, Top2boxes)



This image was created with Generative AI.

Scientifically proven.

Designed to reduce cognitive load.

ZEISS Vision Care collaborated with the University of Tübingen to research whether or not the blur in eyeglass lenses impacts the cognitive load. Researchers surmised that lens blur might make the brain work harder to process visual input.

How did we do that?

We asked participants to perform certain tasks while we monitored their brain activity with EEG (electroencephalography). We then asked them to do these tasks while being impacted by lens blur. Finally, we compared their performance and brain responses in each case.

For the first time worldwide, we demonstrated that lens blur affects the cognitive load.



With this in mind, we designed **ZEISS ClearMind lenses** to reduce the cognitive load and ease the overall mental effort.*

14% less cognitive load compared to a ZEISS standard progressive lens.*

16% higher blink rate on average, suggesting a lower cognitive load compared to a ZEISS standard lens.*

Disclaimer: ZEISS ClearMind lenses do not alleviate the cognitive strain/load caused by stress, fatigue or other psychological factors.
* Performance test on ZEISS ClearMind progressive lenses compared to standard ZEISS progressive lenses, N=32, NASA-TLX assessment of perceived cognitive load in a blinded test. ZEISS Vision Science Lab, Institute for Ophthalmic Research, University of Tuebingen, DE 2025.


ZEISS ClearMind lenses

See clearly. Think clearly.

ZEISS ClearMind lenses incorporate innovative ZEISS NeuroOptix technology, utilizing advanced design intelligence that integrates neuroscience principles to optimize lens design.

How does ZEISS NeuroOptix technology work?

ZEISS NeuroOptix technology aligns the lens clarity zones with natural eye movement patterns, **delivering up to 41% larger clear vision zones.*** This reduces the cognitive load on the brain, providing relaxed, effortless vision and enhanced comfort during long days of digital, visually demanding and dynamic activities.



* than ZEISS standard progressive lenses for intermediate distance. Calculation & Analysis by Technology & Innovation, Carl Zeiss Vision International GmbH, DE 2025 (unpublished, data on file)



How you benefit.

Extremely clear vision.

Enjoy sharper, clearer vision throughout your day thanks to *reduced lens blur*.

96% of wearers confirm extremely clear vision.**

Less cognitive load.

Clear vision helps your brain process visual information more easily, allowing you to *think clearly*.

90% of wearers experience less cognitive load compared to a ZEISS standard lens.***

Better concentration.

Seeing and thinking clearly helps you stay *engaged and focused*, even during demanding tasks.

76% of wearers state that the lenses increased their ability to concentrate on tasks.**

Increased well-being.

Visual clarity make you feel more comfortable, improving *your overall quality of life*.

84% of wearers confirm that the lenses contributed positively to their overall feeling of well-being.**


** Opinion-based survey among wearers of ZEISS ClearMind lenses in Germany, Italy, India & China (N=298, Single Vision, Digital & Progressive lenses), 2025, Carl Zeiss Vision International GmbH, DE (unpublished, data on file, Top2boxes).

*** Performance test on ZEISS ClearMind progressive lenses compared to standard ZEISS progressive lenses, N=32, NASA-TLX assessment of perceived cognitive load in a blinded test. ZEISS Vision Science Lab, Institute for Ophthalmic Research, University of Tuebingen, DE 2025.

ZEISS Single Vision ClearMind lenses

ZEISS Single Vision ClearMind lenses provide extremely clear vision with reduced peripheral blur towards the edges compared to ZEISS standard lenses.

Large zone of clear vision



Optimized for
temporary
blur perception

Low and optimized peripheral blur towards the edges

Up to **48%**
larger zones of
clear vision with
no visible blur
than a ZEISS
standard Single
Vision lenses.¹



Consumer benefits

- Extremely clear vision all day long.
- Vitalized eyes for all-day comfort.
- Reduced mental effort.


1. Calculation & Analysis by Technology & Innovation, Carl Zeiss Vision International GmbH, DE 2025 (unpublished, data on file)

ZEISS Single Vision ClearMind Lenses

See Clearly. Think Clearly.

Extremely clear vision.
Designed to reduce the cognitive load.

ZEISS Single Vision ClearMind lenses are engineered to expand the effective visual field compared to conventional ZEISS single vision lenses by minimizing peripheral blur. This lens design incorporates innovative **ZEISS NeuroOptix technology**, utilizing advanced design intelligence with integrated neuroscience principles. This enhances visual processing efficiency, reducing the cognitive load through uninterrupted visual clarity.




Up to **48%**
larger zones of clear vision with
no visible blur than ZEISS standard
single vision lenses.¹

75% of wearers confirm
that it feels easier for them to stay
focused with these lenses.²



By the way
Vision zones schematically
illustrated here are less
perceptible in reality.



ZEISS Single Vision ClearMind lenses come in 2 tiers:
Superb and Individual 3.
Individual 3 tier considers individual anatomy and special vision needs
to get the most out of the customized lens design for any patient.

¹ Calculation & Analysis by Technology & Innovation, Carl Zeiss Vision International GmbH, DE 2025 (unpublished, data on file)
² Opinion-based survey among wearers of ZEISS ClearMind Lenses in Germany, Italy, India & China (N=298, Single Vision, Digital & Progressive Lenses), 2025, Carl Zeiss Vision International GmbH, DE (unpublished, data on file, Top2boxes)

ZEISS ClearMind Single Vision / Asiana	BlueGuard							
	1.50	1.60	1.67	1.74				
DuraVision® Plus Gold UV	1345	1570	1790	2440				
DuraVision® Plus Platinum UV	1145	1370	1590	2240	Now with BlueGuard			
Clear PhotoFusion X Tinted	UVProtect				PhotoFusion X			
	1.50	1.60	1.67	1.74	1.50	1.60	1.67	1.74**
DuraVision® Plus Mirror UV	1245	1470	-	-	-	-	-	-
DuraVision® Plus Sun UV	1145	1370	1590	-	-	-	-	-
DuraVision® Plus Flash UV	1245	1470	1690	-	1745	1970	2190	-
DuraVision® Plus Gold UV	1245	1470	1690	2340	1845	2070	2290	2940
DuraVision® Plus Platinum UV	1045	1270	1490	2140	1645	1870	2090	2740

* List price for DuraVision® Plus Flash, DuraVision® Plus Sun, DuraVision® Plus Mirror on UV Protect is inclusive of Standard Tint colour, excluding Syklet and Special Tint.

**Top-up RM100/pair for Skylet (Fun, Road and Sport) & Special Tint (ProGolf, ProGolf Gradient and Double Gradient Tints)

**PhotoFusion X 1.74 available in Grey, Extra Grey, Brown and Black only.

AdaptiveSun	AdaptiveSun				AdaptiveSun Polarized		
	1.50	1.60	1.67	1.74	1.50	1.60	1.67
DuraVision® Plus Mirror UV	1945	2170	-	-	2095	2320	-
DuraVision® Plus Sun UV	1845	2070	2290	-	1995	2220	-
DuraVision® Plus Flash UV	1945	2170	2390	-	2095	2320	-
DuraVision® Plus Gold UV	2045	2270	2490	-	2195	2420	-
DuraVision® Plus Platinum UV	1845	2070	2290	-	1995	2220	-
Polarized	1.50	1.60	1.67	1.74			
DuraVision® Plus Mirror UV	1745	1970	-	-			
DuraVision® Plus Sun UV	1645	1870	2090	-			
DuraVision® Plus Flash UV	1745	1970	2190	-			
DuraVision® Plus Gold UV	1845	2070	2290	-			
DuraVision® Plus Platinum UV	1645	1870	2090	-			

* Require customer's date of birth (e.g. 14 Feb 1996), Frame Shape and Size, Pantoscopic angle for lens optimization

* AdaptiveSun Blue is only available in 1.5 and 1.6 index.

* AdaptiveSun and AdaptiveSun Polarized is not available in DuraVision® Plus Mirror UV Strong Blue.

* Tinted colour is required for DuraVision® Plus Mirror UV, DuraVision® Plus Sun UV, DuraVision® Plus Flash UV (except for PhotoFusion X)

BlueGuard / UVProtect	SPH	CYL	COM RX.	DIA	PRISM
1.50	+6.00 to -6.00	-4.00	-6.00	75	3
1.60	+10.00 to -10.00	-6.00	-10.00	75	3
1.67	+10.00 to -17.00	-6.00	-17.00	70	3
1.74	+16.00 to -20.00	-6.00	-20.00	70	3
PhotoFusion X	SPH	CYL	COM RX.	DIA	PRISM
1.50	+6.00 to -6.00	-4.00	-6.00	75	3
1.60	+8.00 to -10.00	-6.00	-10.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	75	3
1.74	+9.00 to -14.00	-6.00	-14.00	75	3
AdaptiveSun	SPH	CYL	COM RX.	DIA	PRISM
1.50	+6.00 to -6.00	-4.00	-6.00	75	3
1.60	+8.00 to -10.00	-6.00	-10.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	70	3
AdaptiveSun Polarized	SPH	CYL	COM RX.	DIA	PRISM
1.50	+6.00 to -6.00	-4.00	-6.00	80	3
1.60	+10.00 to -11.00	-6.00	-11.00	80	3
Polarized	SPH	CYL	COM RX.	DIA	PRISM
1.50	+8.00 to -8.00	-4.00	-8.00	80	3
1.60	+10.00 to -11.00	-6.00	-11.00	80	3
1.67	+11.00 to -12.00	-6.00	-12.00	80	3


* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Price quoted is as per pair in Ringgit Malaysia and subject to prevailing tax. Price is subject to change without prior notice.

ZEISS Digital ClearMind lenses

ZEISS Digital ClearMind lenses provide accommodative support in the lower lens part for comfortable reading of close-up objects.



Consumer benefits

- Extremely clear vision all day long.
- Vitalized eyes that feel less tired after a long day thanks to an extra boost in the near-vision zone that supports close-up tasks.
- Reduced mental effort.

Up to **7%** larger zones of clear vision than ZEISS standard digital lenses for intermediate distance.¹

Up to **28%** larger zones of clear vision than ZEISS standard digital lenses for near distance.¹


Up to **11%** less blur towards the edges compared to ZEISS standard digital lenses¹

ZEISS Digital ClearMind Lenses

See Clearly. Think Clearly.

Extremely clear vision.
Designed to reduce the cognitive load.

ZEISS Digital ClearMind lenses incorporate innovative **ZEISS NeuroOptix technology**, utilizing advanced design intelligence that integrates neuroscience principles to optimize lens design. ZEISS NeuroOptix technology aligns the lens clarity zones with natural eye movement patterns, **delivering large clear vision zones.**¹ This reduces the cognitive load on the brain, providing relaxed, effortless vision and enhanced comfort throughout long days filled with visual stimuli.




By the way
Vision zones schematically illustrated here are less perceptible in reality.



Up to **28%** larger zones of clear vision.²

85% of wearers confirm that the ZEISS Digital ClearMind lenses provide an extra stress reduction for their eyes.³

75% of wearers confirm that it feels easier for them to stay focused with these lenses.



ZEISS Digital ClearMind lenses are designed for people experiencing early signs of aging vision and are good practice for eventually transitioning to progressive lenses.

ZEISS Digital ClearMind lenses come in 2 tiers: Superb and Individual 3. Individual 3 tier considers individual anatomy and special vision needs to get the most out of the customized lens design for any patient.

1 Calculation & Analysis by Technology & Innovation, Carl Zeiss Vision International GmbH, DE 2025 (unpublished, data on file)
2 Compared to ZEISS standard digital lenses for near distance. Calculation & Analysis by Technology & Innovation, Carl Zeiss Vision International GmbH, DE 2025 (unpublished, data on file)
3 Opinion-based survey among wearers of ZEISS ClearMind lenses in Germany, Italy, India & China (N=64, Digital), 2025, Carl Zeiss Vision International GmbH, DE (unpublished, data on file, Top2boxes)
4 Opinion-based survey among wearers of ZEISS ClearMind lenses in Germany, Italy, India & China (N=298, Single Vision, Digital & Progressive lenses), 2025, Carl Zeiss Vision International GmbH, DE (unpublished, data on file, Top2boxes)

ZEISS ClearMind Digital / Asiana	BlueGuard				Now with BlueGuard			
	1.50	1.60	1.67	1.74				
DuraVision® Plus Gold UV	1455	1790	2140	2950				
DuraVision® Plus Platinum UV	1255	1590	1940	2750				
Clear PhotoFusion X Tinted	UVProtect				PhotoFusion X			
	1.50	1.60	1.67	1.74	1.50	1.60	1.67	1.74**
DuraVision® Plus Mirror UV	1355	1690	-	-	-	-	-	-
DuraVision® Plus Sun UV	1255	1590	1940	-	-	-	-	-
DuraVision® Plus Flash UV	1355	1690	2040	-	1855	2190	2540	-
DuraVision® Plus Gold UV	1355	1690	2040	2850	1955	2290	2640	3450
DuraVision® Plus Platinum UV	1155	1490	1840	2650	1755	2090	2440	3250

* List price for DuraVision® Plus Flash, DuraVision® Plus Sun, DuraVision® Plus Mirror on UV Protect is inclusive of Standard Tint colour, excluding Syklet and Special Tint.

** Top-up RM100/pair for Skylet (Fun, Road and Sport) & Special Tint (ProGolf, ProGolf Gradient and Double Gradient Tints)

** PhotoFusion X 1.74 available in Grey, Extra Grey, Brown and Black only.

AdaptiveSun	AdaptiveSun				AdaptiveSun Polarized		
	1.50	1.60	1.67	1.74	1.50	1.60	1.67
DuraVision® Plus Mirror UV	2055	2390	-	-	2205	2540	-
DuraVision® Plus Sun UV	1955	2290	2640	-	2105	2440	-
DuraVision® Plus Flash UV	2055	2390	2740	-	2205	2540	-
DuraVision® Plus Gold UV	2155	2490	2840	-	2305	2640	-
DuraVision® Plus Platinum UV	1955	2290	2640	-	2105	2440	-
Polarized	1.50	1.60	1.67	1.74			
DuraVision® Plus Mirror UV	1855	2190	-	-			
DuraVision® Plus Sun UV	1755	2090	2440	-			
DuraVision® Plus Flash UV	1855	2190	2540	-			
DuraVision® Plus Gold UV	1955	2290	2640	-			
DuraVision® Plus Platinum UV	1755	2090	2440	-			

* Require customer's date of birth (e.g. 14 Feb 1996). Frame Shape and Size for lens optimization

* AdaptiveSun Blue is only available in 1.5 and 1.6 index.

* AdaptiveSun and AdaptiveSun Polarized is not available in DuraVision® Plus Mirror UV Strong Blue.

* Tinted colour is required for DuraVision® Plus Mirror UV, DuraVision® Plus Sun UV, DuraVision® Plus Flash UV (except for PhotoFusion)

BlueGuard / UVProtect	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+6.00 to -7.00	-6.00	-7.00	+0.50 to +1.25	75	3
1.60	+10.00 to -10.00	-6.00	-10.00	+0.50 to +1.25	75	3
1.67	+8.00 to -17.00	-6.00	-17.00	+0.50 to +1.25	70	3
1.74	+16.00 to -20.00	-6.00	-20.00	+0.50 to +1.25	70	3
PhotoFusion X	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.50 to +1.25	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.50 to +1.25	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.50 to +1.25	75	3
1.74	+9.00 to -14.00	-6.00	-14.00	+0.50 to +1.25	75	3
AdaptiveSun	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.50 to +1.25	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.50 to +1.25	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.50 to +1.25	70	3
AdaptiveSun Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.50 to +1.25	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.50 to +1.25	80	3
Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.50 to +1.25	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.50 to +1.25	80	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.50 to +1.25	80	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

* You are advised to use VISUSTORE® to verify the availability of maximum diameter.


Price quoted is as per pair in Ringgit Malaysia and subject to prevailing tax. Price is subject to change without prior notice.

ZEISS Progressive ClearMind lenses

Clarity without compromise:
51% adapted in an hour.
88% within days.²

ZEISS Progressive ClearMind lenses have a hybrid optimized design that provides large fields of clear vision at all distances and in all directions

This design features an extra large intermediate zone and very low, optimized peripheral blur. These elements support a smooth transition between blur zones, allowing for fast adaptation and extremely clear vision at any distance and direction for the wearer.



Consumer benefits

- Extremely clear vision all day long.
- Vitalized eyes for all day comfort
- Smooth transitions between the near and far vision zones of the lens.
- Fast adaptation.
- Reduced mental effort.

Up to **10%** less blur towards the edges compared to ZEISS SmartLife progressive lenses

Up to **9%** larger zones of clear vision than ZEISS SmartLife progressive lenses for far distance.

Up to **26%** larger zones of clear vision than ZEISS SmartLife progressive lenses for intermediate distance.

Up to **41%** larger zones of clear vision than ZEISS SmartLife progressive lenses for near distance.

1. Calculation & Analysis by Technology & Innovation, Carl Zeiss Vision International GmbH, DE 2025 (unpublished, data on file).


2. Opinion-based survey among wearers of ZEISS ClearMind lenses in Germany, Italy, India & China (N=193 Progressive lenses), 2025, Carl Zeiss Vision International GmbH, DE (unpublished, data on file, Top2boxes)

ZEISS Progressive ClearMind Lenses

See Clearly. Think Clearly.

Extremely clear vision.
Designed to reduce the cognitive load.

ZEISS Progressive ClearMind lenses incorporate innovative **ZEISS NeurOptix technology**, utilizing advanced design intelligence that integrates neuroscience principles to optimize lens design. ZEISS NeurOptix technology aligns the lens clarity zones with natural eye movement patterns, **delivering large clear vision zones**.¹ This reduces the cognitive load on the brain, providing relaxed, effortless vision and enhanced comfort throughout long days filled with visual stimuli.



ZEISS ClearMind Progressive Superb	BlueGuard							
	1.50	1.60	1.67	1.74				
DuraVision® Plus Gold UV	3950	4480	5030	5580				
DuraVision® Plus Platinum UV	3750	4280	4830	5380	Now with BlueGuard			
Clear PhotoFusion X Tinted	UVProtect				PhotoFusion X			
	1.50	1.60	1.67	1.74	1.50	1.60	1.67	1.74**
DuraVision® Plus Mirror UV	3850	4380	-	-	-	-	-	-
DuraVision® Plus Sun UV	3750	4280	4830	-	-	-	-	-
DuraVision® Plus Flash UV	3850	4380	4930	-	4350	4880	5430	-
DuraVision® Plus Gold UV	3850	4380	4930	5480	4450	4980	5530	6080
DuraVision® Plus Platinum UV	3650	4180	4730	5280	4250	4780	5330	5880

* List price for DuraVision® Plus Flash, DuraVision® Plus Sun, DuraVision® Plus Mirror on UV Protect is inclusive of Standard Tint colour, excluding Syklet and Special Tint.

* Top-up RM100/pair for Skylet (Fun, Road and Sport) & Special Tint (ProGolf, ProGolf Gradient and Double Gradient Tints)

**PhotoFusion X 1.74 available in Grey, Extra Grey, Brown and Black only.

AdaptiveSun	AdaptiveSun				AdaptiveSun Polarized		
	1.50	1.60	1.67	1.74	1.50	1.60	1.67
DuraVision® Plus Mirror UV	4550	5080	-	-	4700	5230	-
DuraVision® Plus Sun UV	4450	4980	5530	-	4600	5130	-
DuraVision® Plus Flash UV	4550	5080	5630	-	4700	5230	-
DuraVision® Plus Gold UV	4650	5180	5730	-	4800	5330	-
DuraVision® Plus Platinum UV	4450	4980	5530	-	4600	5130	-
Polarized	1.50	1.60	1.67	1.74			
DuraVision® Plus Mirror UV	4350	4880	-	-			
DuraVision® Plus Sun UV	4250	4780	5330	-			
DuraVision® Plus Flash UV	4350	4880	5430	-			
DuraVision® Plus Gold UV	4450	4980	5530	-			
DuraVision® Plus Platinum UV	4250	4780	5330	-			

* Require customer's date of birth (e.g. 14 Feb 1970), Frame Shape and Size, Pantoscopic angle, Vertex Distance, FrameFit value (FFV from -1.0 to 6.0) for lens optimization

* AdaptiveSun Blue is only available in 1.5 and 1.6 index.

* AdaptiveSun and AdaptiveSun Polarized is not available in DuraVision® Plus Mirror UV Strong Blue.

* Tinted colour is required for DuraVision® Plus Mirror UV, DuraVision® Plus Sun UV, DuraVision® Plus Flash UV (except for PhotoFusion X)

BlueGuard / UVProtect	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+6.00 to -7.00	-6.00	-7.00	+0.75 to +3.50	75	3
1.60	+10.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -17.00	-6.00	-17.00	+0.75 to +3.50	70	3
1.74	+16.00 to -20.00	-6.00	-20.00	+0.75 to +3.50	70	3
PhotoFusion X	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	75	3
1.74	+9.00 to -14.00	-6.00	-14.00	+0.75 to +3.50	75	3
AdaptiveSun	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	70	3
AdaptiveSun Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.75 to +3.50	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.75 to +3.50	80	3
Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.75 to +3.50	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.75 to +3.50	80	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	80	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Price quoted is as per pair in Ringgit Malaysia and subject to prevailing tax. Price is subject to change without prior notice.

ZEISS ClearMind Progressive Plus	BlueGuard				Now with BlueGuard			
	1.50	1.60	1.67	1.74				
DuraVision® Plus Gold UV	2950	3440	3990	4480				
DuraVision® Plus Platinum UV	2750	3240	3790	4280				
Clear PhotoFusion X Tinted	UVProtect				PhotoFusion X			
	1.50	1.60	1.67	1.74	1.50	1.60	1.67	1.74**
DuraVision® Plus Mirror UV	2850	3340	-	-	-	-	-	-
DuraVision® Plus Sun UV	2750	3240	3790	-	-	-	-	-
DuraVision® Plus Flash UV	2850	3340	3890	-	3350	3840	4390	-
DuraVision® Plus Gold UV	2850	3340	3890	4380	3450	3940	4490	4980
DuraVision® Plus Platinum UV	2650	3140	3690	4180	3250	3740	4290	4780

* List price for DuraVision® Plus Flash, DuraVision® Plus Sun, DuraVision® Plus Mirror on UV Protect is inclusive of Standard Tint colour, excluding Skyket and Special Tint.

* Top-up RM100/pair for Skyket (Fun, Road and Sport) & Special Tint (ProGolf, ProGolf Gradient and Double Gradient Tints)

**PhotoFusion X 1.74 available in Grey, Extra Grey, Brown and Black only.

AdaptiveSun	AdaptiveSun				AdaptiveSun Polarized		
	1.50	1.60	1.67	1.74	1.50	1.60	1.67
DuraVision® Plus Mirror UV	3550	4040	-	-	3700	4190	-
DuraVision® Plus Sun UV	3450	3940	4490	-	3600	4090	-
DuraVision® Plus Flash UV	3550	4040	4590	-	3700	4190	-
DuraVision® Plus Gold UV	3650	4140	4690	-	3800	4290	-
DuraVision® Plus Platinum UV	3450	3940	4490	-	3600	4090	-
Polarized	1.50	1.60	1.67	1.74			
DuraVision® Plus Mirror UV	3350	3840	-	-			
DuraVision® Plus Sun UV	3250	3740	4290	-			
DuraVision® Plus Flash UV	3350	3840	4390	-			
DuraVision® Plus Gold UV	3450	3940	4490	-			
DuraVision® Plus Platinum UV	3250	3740	4290	-			

* Require customer's date of birth (e.g. 14 Feb 1970), Frame Shape and Size, FrameFit value (FFV from -1.0 to 6.0) for lens optimization

* AdaptiveSun Blue is only available in 1.5 and 1.6 index.

* AdaptiveSun and AdaptiveSun Polarized is not available in DuraVision® Plus Mirror UV Strong Blue.

* Tinted colour is required for DuraVision® Plus Mirror UV, DuraVision® Plus Sun UV, DuraVision® Plus Flash UV (except for PhotoFusion X)

BlueGuard / UVProtect	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+6.00 to -7.00	-6.00	-7.00	+0.75 to +3.50	75	3
1.60	+10.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -17.00	-6.00	-17.00	+0.75 to +3.50	70	3
1.74	+16.00 to -20.00	-6.00	-20.00	+0.75 to +3.50	70	3
PhotoFusion X	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	75	3
1.74	+9.00 to -14.00	-6.00	-14.00	+0.75 to +3.50	75	3
AdaptiveSun	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	70	3
AdaptiveSun Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.75 to +3.50	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.75 to +3.50	80	3
Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.75 to +3.50	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.75 to +3.50	80	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	80	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Price quoted is as per pair in Ringgit Malaysia and subject to prevailing tax. Price is subject to change without prior notice.

ZEISS ClearMind Progressive Pure	BlueGuard							
	1.50	1.60	1.67	1.74				
DuraVision® Plus Gold UV	2350	2750	3170	3800				
DuraVision® Plus Platinum UV	2150	2550	2970	3600	Now with BlueGuard			
Clear PhotoFusion X Tinted	UVProtect				PhotoFusion X			
	1.50	1.60	1.67	1.74	1.50	1.60	1.67	1.74**
DuraVision® Plus Mirror UV	2250	2650	-	-	-	-	-	-
DuraVision® Plus Sun UV	2150	2550	2970	-	-	-	-	-
DuraVision® Plus Flash UV	2250	2650	3070	-	2750	3150	3570	-
DuraVision® Plus Gold UV	2250	2650	3070	3700	2850	3250	3670	4300
DuraVision® Plus Platinum UV	2050	2450	2870	3500	2650	3050	3470	4100

* List price for DuraVision® Plus Flash, DuraVision® Plus Sun, DuraVision® Plus Mirror on UV Protect is inclusive of Standard Tint colour, excluding Syklet and Special Tint.

* Top-up RM100/pair for Skylet (Fun, Road and Sport) & Special Tint (ProGolf, ProGolf Gradient and Double Gradient Tints)

**PhotoFusion X 1.74 available in Grey, Extra Grey, Brown and Black only.

AdaptiveSun	AdaptiveSun				AdaptiveSun Polarized		
	1.50	1.60	1.67	1.74	1.50	1.60	1.67
DuraVision® Plus Mirror UV	2950	3350	-	-	3100	3500	-
DuraVision® Plus Sun UV	2850	3250	3670	-	3000	3400	-
DuraVision® Plus Flash UV	2950	3350	3770	-	3100	3500	-
DuraVision® Plus Gold UV	3050	3450	3870	-	3200	3600	-
DuraVision® Plus Platinum UV	2850	3250	3670	-	3000	3400	-
Polarized	1.50	1.60	1.67	1.74			
DuraVision® Plus Mirror UV	2750	3150	-	-			
DuraVision® Plus Sun UV	2650	3050	3470	-			
DuraVision® Plus Flash UV	2750	3150	3570	-			
DuraVision® Plus Gold UV	2850	3250	3670	-			
DuraVision® Plus Platinum UV	2650	3050	3470	-			

* Require customer's date of birth (e.g. 14 Feb 1970), Frame Shape and Size for lens optimization

* Require specific lens corridor: Regular (14mm), Medium (12mm) or Short (10mm).

* AdaptiveSun Blue is only available in 1.5 and 1.6 index.

* AdaptiveSun and AdaptiveSun Polarized is not available in DuraVision® Plus Mirror UV Strong Blue.

* Tinted colour is required for DuraVision® Plus Mirror UV, DuraVision® Plus Sun UV, DuraVision® Plus Flash UV (except for PhotoFusion X)

BlueGuard / UVProtect	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+6.00 to -7.00	-6.00	-7.00	+0.75 to +3.50	75	3
1.60	+10.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -17.00	-6.00	-17.00	+0.75 to +3.50	70	3
1.74	+16.00 to -20.00	-6.00	-20.00	+0.75 to +3.50	70	3
PhotoFusion X	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	75	3
1.74	+9.00 to -14.00	-6.00	-14.00	+0.75 to +3.50	75	3
AdaptiveSun	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	70	3
AdaptiveSun Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.75 to +3.50	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.75 to +3.50	80	3
Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.75 to +3.50	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.75 to +3.50	80	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	80	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Price quoted is as per pair in Ringgit Malaysia and subject to prevailing tax. Price is subject to change without prior notice.



ZEISS Mineral Lenses | The material you choose exceeds your expectations.


The unique structure of the mineral material offers you the best protection in our daily work life and prevents fine scratches in dusty environments.

The most important qualities of mineral lens material :

- A distinct surface hardness
- An unrivaled scratch resistance leading to an excellent longevity
- Extreme thin lens materials compared to plastic lenses
- More clarity due to higher Abbe numbers and great optical performance

Impressive scratch resistance


Mineral materials guarantee a long lasting life due to their highly crosslinked and dense 3D structure. This is a benefit for daily work life and prevents fine scratches in dusty environments. A clearer vision by looking through the lenses lets your patients enjoy every moment. The mineral material will truly fulfill your patients' need for a long and satisfied usage of their glasses.



Scheme of structured mineral glass molecules 1.6 (left) compared to unstructured plastic molecules (right).

Extreme thin lens materials

Materials with indices up to 1.8 guarantee a thin surface edge for lenses. This ensures a more aesthetic appearance for all consumers that require a thin lens solution. A major benefit especially for consumers with natural high prescriptions who experience presbyopia. A thinner material also guarantees a reduced appearance of the downsizing effect of the eyes with high minus power lenses. All ZEISS mineral progressive lenses are manufactured with the latest freeform technology for a more accurate surfacing than ever before.



Mineral 1.6 (left) compared to mineral 1.8 (right), -7.0 D and @ 70 mm of illustrated lenses.



Scan here for more details
ZEISS Mineral Lenses



Seeing beyond

ZEISS ClearMind Progressive Superb Mineral Lens	Clear		Umbramatic
	1.60	1.80	1.60
Super ET	4280	5380	4880

* Require customer's date of birth (e.g. 14 Feb 1970), Frame Shape and Size, Pantoscopic Angle, Back Vertex Distance, FrameFit (FFV: 0.0 - 6.0) for lens optimization.

* Require 14 days delivery.

* Tinted colour is unavailable.

Clear	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.60	+6.50 to -11.00	-6.00	-11.00	+0.75 to +3.50	70	3
1.80	0.00 to -20.00	-8.00	-20.00	+0.75 to +3.50	65	3
Umbramatic	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.60	+6.50 to -9.00	-4.00	-9.00	+0.75 to +3.50	70	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

You are advised to use VISUSTORE® to verify the availability of maximum diameter.

ZEISS ClearMind Progressive Plus	Clear		Umbramatic
	1.60	1.80	1.60
Super ET	3240	4280	3840

* Require customer's date of birth (e.g. 14 Feb 1970), Frame Shape and Size, Pantoscopic Angle, Back Vertex Distance, FrameFit (FFV: 0.0 - 6.0) for lens optimization.

* Require 14 days delivery.

* Tinted colour is unavailable.

Clear	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.60	+6.50 to -11.00	-6.00	-11.00	+0.75 to +3.50	70	3
1.80	0.00 to -20.00	-8.00	-20.00	+0.75 to +3.50	65	3
Umbramatic	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.60	+6.50 to -9.00	-4.00	-9.00	+0.75 to +3.50	70	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

You are advised to use VISUSTORE® to verify the availability of maximum diameter.

ZEISS ClearMind Digital Lens	Clear		Umbramatic
	1.60	1.80	1.60
Super ET	1590	2750	2190

* Require customer's date of birth (e.g. 14 Feb 1970), Frame Shape and Size, Pantoscopic Angle, Back Vertex Distance, FrameFit (FFV: 0.0 - 6.0) for lens optimization.

* Require 14 days delivery.


* Tinted colour is unavailable.

Clear	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.60	+6.50 to -11.00	-6.00	-11.00	+0.50 to +1.25	70	3
1.80	0.00 to -20.00	-8.00	-20.00	+0.50 to +1.25	65	3
Umbramatic	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.60	+6.50 to -9.00	-4.00	-9.00	+0.50 to +1.25	70	3


* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Stamp and Engravings




ZEISS ClearMind Single Vision lenses




ZEISS Digital ClearMind Lenses lenses

Stamp and Engravings



ZEISS Progressive ClearMind Plus / Superb lens



ZEISS Progressive ClearMind Pure lens

ZEISS ClearMind Lenses

Lens Engraving

Nasal Engraving	ClearMind
ZEISS Single Vision ClearMind	VCMxx
ZEISS Digital ClearMind	BCMxx
ZEISS Digital Asiana* ClearMind	DCAxx
ZEISS Progressive ClearMind Pure	CPUxx
ZEISS Progressive ClearMind Plus	CMPxx
ZEISS Progressive ClearMind Plus Asiana*	CPAxx
ZEISS Progressive ClearMind Superb	CMSxx

*Asiana version only available for some countries.