

Seeing beyond

ZEISS

# ZEISS ClearMind Individual 3

Experience truly personalised vision correction with the Individual 3 customisation level — crafted exclusively for your unique perspective.

# See the world your way.

## Lenses customized especially for your eyes.

Unlock the full potential of your vision with the **Individual 3 customization level** for ZEISS ClearMind lenses. We use ZEISS Intelligence Augmented Design technology to enhance your personal, anatomical and physiological data with **12.5 million data points** about real-life, age-related visual behaviour collected from by ZEISS Global Vision Study.\*

Thoughtfully designed to suit your unique anatomy, specific vision needs and age-related visual behaviour, these lenses adapt perfectly to your favourite frames to deliver extremely clear vision and comfort.

Our partner eye care professionals use the latest ZEISS measuring technology to personalise the perfect ZEISS ClearMind lenses for you and your eyes.




\* ZEISS Global Vision Study 2020/21.

\*\* for intermediate distance (with ZEISS Progressive ClearMind Individual 3 compared to ZEISS Progressive ClearMind Superb Lenses). Calculation & Analysis by Technology & Innovation, Carl Zeiss Vision International GmbH, DE 2025 (unpublished, data on file)

# ZEISS Digital & Progressive ClearMind lenses

Tailored solutions to serve specific needs

| DIGITAL     | 35+ | For individuals 35+ who experience blurry vision, near or far, along with increased eye fatigue due to age-related changes. <b>Accommodative support ranging from +0.50D to 1.25D.</b> | Lens options                                  | Vision zones                                                                                                                                                                                                                                                                                                                                                                                                   | Frame choice                                                                                                                                                                                             |
|-------------|-----|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|             |     |                                                                                                                                                                                        | OUR <b>BEST</b> SOLUTION<br>Individual 3 tier | <ul style="list-style-type: none"><li>• Larger clear zones when fitted to individual anatomy</li><li>• Accommodative support in the lower lens part for relaxed reading Smooth peripheral blur zones</li><li>• Optimised for digital device use</li><li>• Optimised for individuals wearing addition power lenses for the first time</li><li>• Optimised for support of the cognitive load reduction</li></ul> | <br><b>For any frame</b> of your choice, adapted to <b>your individual anatomy</b> and <b>special vision needs.</b>   |
|             |     |                                                                                                                                                                                        | OUR <b>BETTER</b> SOLUTION<br>Superb tier     | <ul style="list-style-type: none"><li>• Large clear zones</li><li>• Accommodative support in the lower lens part for relaxed reading Smooth peripheral blur zones</li><li>• Optimised for digital device use</li><li>• Optimized for individuals wearing addition power lenses for the first time</li><li>• Optimised for support of the cognitive load reduction</li></ul>                                    | <br><b>For any frame</b> of your choice, adapted to <b>your individual anatomy.</b>                                   |
|             |     |                                                                                                                                                                                        |                                               |                                                                                                                                                                                                                                                                                                                                                                                                                |                                                                                                                                                                                                          |
| PROGRESSIVE | 45+ | For individuals 45+ who need seamless vision from near to far and at all distances. <b>From +0.75D up to +4.00D addition power.</b>                                                    | OUR <b>BEST</b> SOLUTION<br>Individual 3 tier | <ul style="list-style-type: none"><li>• Larger clear zones than ZEISS</li><li>• Progressive ClearMind Superb when fitted to individual anatomy &amp; vision needs. Smooth peripheral blur zones.</li><li>• Optimised for digital device use.</li><li>• Optimised for support of the cognitive load reduction.</li></ul> <div>Up to <b>19%</b> larger vision zones.*</div>                                      | <br><b>For any frame</b> of your choice, adapted to <b>your individual anatomy</b> and <b>special vision needs.</b> |
|             |     |                                                                                                                                                                                        | OUR <b>BETTER</b> SOLUTION<br>Superb tier     | <ul style="list-style-type: none"><li>• Larger clear zones than ZEISS</li><li>• Progressive ClearMind Plus when fitted to individual anatomy.</li><li>• Smooth peripheral blur.</li><li>• Optimised for digital device use.</li><li>• Optimised for support of the cognitive load reduction.</li></ul>                                                                                                         | <br><b>For any frame</b> of your choice, adapted to <b>your individual anatomy.</b>                                 |
|             |     |                                                                                                                                                                                        | OUR <b>GOOD</b> SOLUTION<br>Plus tier         | <ul style="list-style-type: none"><li>• Larger clear zones than ZEISS</li><li>• Progressive ClearMind Pure when fitted in near vision zone.</li><li>• Smooth peripheral blur.</li><li>• Optimised for digital device use.</li><li>• Optimised for support of the cognitive load reduction.</li></ul>                                                                                                           | <br><b>For any frame</b> with variable fitting heights up from 13 mm.                                               |
|             |     |                                                                                                                                                                                        | OUR <b>BASIC</b> SOLUTION<br>Pure tier        | <ul style="list-style-type: none"><li>• Large clear zones.</li><li>• Smooth peripheral blur.</li><li>• Optimised for digital device usage.</li><li>• Optimised for support of the cognitive load reduction.</li></ul>                                                                                                                                                                                          | <br><b>For standard frames</b> available in 3 fitting heights (14/16/18 mm).                                        |
|             |     |                                                                                                                                                                                        |                                               |                                                                                                                                                                                                                                                                                                                                                                                                                |                                                                                                                                                                                                          |

\* For intermediate distance (with ZEISS Progressive ClearMind Individual 3 compared to ZEISS Progressive ClearMind Superb Lenses). Calculation & Analysis by Technology & Innovation, Carl Zeiss Vision International GmbH, DE 2025 (unpublished, data on file)

| ZEISS ClearMind<br>Single Vision Individual 3 | BlueGuard |      |      |      |                    |      |      |        |
|-----------------------------------------------|-----------|------|------|------|--------------------|------|------|--------|
|                                               | 1.50      | 1.60 | 1.67 | 1.74 |                    |      |      |        |
| DuraVision® Plus Gold UV                      | 1980      | 2230 | 2470 | 3250 |                    |      |      |        |
| DuraVision® Plus Platinum UV                  | 1780      | 2030 | 2270 | 3050 | Now with BlueGuard |      |      |        |
| Clear   PhotoFusion X   Tinted                | UVProtect |      |      |      | PhotoFusion X      |      |      |        |
|                                               | 1.50      | 1.60 | 1.67 | 1.74 | 1.50               | 1.60 | 1.67 | 1.74** |
| DuraVision® Plus Mirror UV                    | 1880      | 2130 | -    | -    | -                  | -    | -    | -      |
| DuraVision® Plus Sun UV                       | 1780      | 2030 | 2270 | -    | -                  | -    | -    | -      |
| DuraVision® Plus Flash UV                     | 1880      | 2130 | 2370 | -    | 2380               | 2630 | 2870 | -      |
| DuraVision® Plus Gold UV                      | 1880      | 2130 | 2370 | 3150 | 2480               | 2730 | 2970 | 3750   |
| DuraVision® Plus Platinum UV                  | 1680      | 1930 | 2170 | 2950 | 2280               | 2530 | 2770 | 3550   |

\* List price for DuraVision® Plus Flash, DuraVision® Plus Sun, DuraVision® Plus Mirror on UV Protect is inclusive of Standard Tint colour, excluding Skyket and Special Tint.

\* Top-up RM100/pair for Skyket (Fun, Road and Sport) & Special Tint (ProGolf, ProGolf Gradient and Double Gradient Tints)

\*\*PhotoFusion X 1.74 available in Grey, Extra Grey, Brown and Black only.

| AdaptiveSun                  | AdaptiveSun |      |      |      | AdaptiveSun Polarized |      |      |
|------------------------------|-------------|------|------|------|-----------------------|------|------|
|                              | 1.50        | 1.60 | 1.67 | 1.74 | 1.50                  | 1.60 | 1.67 |
| DuraVision® Plus Mirror UV   | 2580        | 2830 | -    | -    | 2730                  | 2980 | -    |
| DuraVision® Plus Sun UV      | 2480        | 2730 | 2970 | -    | 2630                  | 2880 | -    |
| DuraVision® Plus Flash UV    | 2580        | 2830 | 3070 | -    | 2730                  | 2980 | -    |
| DuraVision® Plus Gold UV     | 2680        | 2930 | 3170 | -    | 2830                  | 3080 | -    |
| DuraVision® Plus Platinum UV | 2480        | 2730 | 2970 | -    | 2630                  | 2880 | -    |
| Polarized                    | 1.50        | 1.60 | 1.67 | 1.74 |                       |      |      |
| DuraVision® Plus Mirror UV   | 2380        | 2630 | -    | -    |                       |      |      |
| DuraVision® Plus Sun UV      | 2280        | 2530 | 2770 | -    |                       |      |      |
| DuraVision® Plus Flash UV    | 2380        | 2630 | 2870 | -    |                       |      |      |
| DuraVision® Plus Gold UV     | 2480        | 2730 | 2970 | -    |                       |      |      |
| DuraVision® Plus Platinum UV | 2280        | 2530 | 2770 | -    |                       |      |      |

\* Require customer's date of birth (e.g. 14 Feb 1996). Frame Shape and Size, Pantoscopic angle, Wrap angle and Back Vertex Distance, for lens optimization

\* AdaptiveSun Blue is only available in 1.5 and 1.6 index.

\* AdaptiveSun and AdaptiveSun Polarized is not available in DuraVision® Plus Mirror UV Strong Blue.

\* Tinted colour is required for DuraVision® Plus Mirror UV, DuraVision® Plus Sun UV, DuraVision® Plus Flash UV (except for PhotoFusion X)

| BlueGuard / UVProtect | SPH              | CYL   | COM RX. | DIA | PRISM |
|-----------------------|------------------|-------|---------|-----|-------|
| 1.50                  | +6.00 to -6.00   | -4.00 | -6.00   | 75  | 3     |
| 1.60                  | +10.00 to -10.00 | -6.00 | -10.00  | 75  | 3     |
| 1.67                  | +10.00 to -17.00 | -6.00 | -17.00  | 70  | 3     |
| 1.74                  | +16.00 to -20.00 | -6.00 | -20.00  | 70  | 3     |
| PhotoFusion X         | SPH              | CYL   | COM RX. | DIA | PRISM |
| 1.50                  | +6.00 to -6.00   | -4.00 | -6.00   | 75  | 3     |
| 1.60                  | +8.00 to -10.00  | -6.00 | -10.00  | 75  | 3     |
| 1.67                  | +8.00 to -12.00  | -6.00 | -12.00  | 70  | 3     |
| 1.74                  | +9.00 to -14.00  | -6.00 | -14.00  | 75  | 3     |
| AdaptiveSun           | SPH              | CYL   | COM RX. | DIA | PRISM |
| 1.50                  | +6.00 to -6.00   | -4.00 | -6.00   | 75  | 3     |
| 1.60                  | +8.00 to -10.00  | -6.00 | -10.00  | 75  | 3     |
| 1.67                  | +8.00 to -12.00  | -6.00 | -12.00  | 70  | 3     |
| AdaptiveSun Polarized | SPH              | CYL   | COM RX. | DIA | PRISM |
| 1.50                  | +6.00 to -6.00   | -4.00 | -6.00   | 80  | 3     |
| 1.60                  | +10.00 to -11.00 | -6.00 | -11.00  | 80  | 3     |
| Polarized             | SPH              | CYL   | COM RX. | DIA | PRISM |
| 1.50                  | +8.00 to -8.00   | -4.00 | -8.00   | 80  | 3     |
| 1.60                  | +10.00 to -11.00 | -6.00 | -11.00  | 80  | 3     |
| 1.67                  | +11.00 to -12.00 | -6.00 | -12.00  | 80  | 3     |

\* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

\* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Price quoted is as per pair in Ringgit Malaysia and subject to prevailing tax. Price is subject to change without prior notice.

| ZEISS ClearMind<br>Digital Individual 3 | BlueGuard |      |      |      |                    |      |      |        |
|-----------------------------------------|-----------|------|------|------|--------------------|------|------|--------|
|                                         | 1.50      | 1.60 | 1.67 | 1.74 |                    |      |      |        |
| DuraVision® Plus Gold UV                | 2150      | 2390 | 2800 | 3270 |                    |      |      |        |
| DuraVision® Plus Platinum UV            | 1950      | 2190 | 2600 | 3070 | Now with BlueGuard |      |      |        |
| Clear   PhotoFusion X   Tinted          | UVProtect |      |      |      | PhotoFusion X      |      |      |        |
|                                         | 1.50      | 1.60 | 1.67 | 1.74 | 1.50               | 1.60 | 1.67 | 1.74** |
| DuraVision® Plus Mirror UV              | 2050      | 2290 | -    | -    | -                  | -    | -    | -      |
| DuraVision® Plus Sun UV                 | 1950      | 2190 | 2600 | -    | -                  | -    | -    | -      |
| DuraVision® Plus Flash UV               | 2050      | 2290 | 2700 | -    | 2550               | 2790 | 3200 | -      |
| DuraVision® Plus Gold UV                | 2050      | 2290 | 2700 | 3170 | 2650               | 2890 | 3300 | 3770   |
| DuraVision® Plus Platinum UV            | 1850      | 2090 | 2500 | 2970 | 2450               | 2690 | 3100 | 3570   |

\* List price for DuraVision® Plus Flash, DuraVision® Plus Sun, DuraVision® Plus Mirror on UV Protect is inclusive of Standard Tint colour, excluding Syklet and Special Tint.

\* Top-up RM100/pair for Syklet (Fun, Road and Sport) & Special Tint (ProGolf, ProGolf Gradient and Double Gradient Tints)

\*\* PhotoFusion X 1.74 available in Grey, Extra Grey, Brown and Black only.

| AdaptiveSun                  | AdaptiveSun |      |      |      | AdaptiveSun Polarized |      |      |
|------------------------------|-------------|------|------|------|-----------------------|------|------|
|                              | 1.50        | 1.60 | 1.67 | 1.74 | 1.50                  | 1.60 | 1.67 |
| DuraVision® Plus Mirror UV   | 2750        | 2990 | -    | -    | 2900                  | 3140 | -    |
| DuraVision® Plus Sun UV      | 2650        | 2890 | 3300 | -    | 2800                  | 3040 | -    |
| DuraVision® Plus Flash UV    | 2750        | 2990 | 3400 | -    | 2900                  | 3140 | -    |
| DuraVision® Plus Gold UV     | 2850        | 3090 | 3500 | -    | 3000                  | 3240 | -    |
| DuraVision® Plus Platinum UV | 2650        | 2890 | 3300 | -    | 2800                  | 3040 | -    |
| Polarized                    | 1.50        | 1.60 | 1.67 | 1.74 |                       |      |      |
| DuraVision® Plus Mirror UV   | 2550        | 2790 | -    | -    |                       |      |      |
| DuraVision® Plus Sun UV      | 2450        | 2690 | 3100 | -    |                       |      |      |
| DuraVision® Plus Flash UV    | 2550        | 2790 | 3200 | -    |                       |      |      |
| DuraVision® Plus Gold UV     | 2650        | 2890 | 3300 | -    |                       |      |      |
| DuraVision® Plus Platinum UV | 2450        | 2690 | 3100 | -    |                       |      |      |

\* Require customer's date of birth (eg. 14 Feb 1996), Frame Shape and Size, Pantoscopic angle for lens optimization, FFV

\* AdaptiveSun Blue is only available in 1.5 and 1.6 index.

\* AdaptiveSun and AdaptiveSun Polarized is not available in DuraVision® Plus Mirror UV Strong Blue.

\* Tinted colour is required for DuraVision® Plus Mirror UV, DuraVision® Plus Sun UV, DuraVision® Plus Flash UV (except for PhotoFusion X)

| BlueGuard / UVProtect | SPH              | CYL   | COM RX. | ADD            | DIA | PRISM |
|-----------------------|------------------|-------|---------|----------------|-----|-------|
| 1.50                  | +6.00 to -7.00   | -6.00 | -7.00   | +0.50 to +1.25 | 75  | 3     |
| 1.60                  | +10.00 to -10.00 | -6.00 | -10.00  | +0.50 to +1.25 | 75  | 3     |
| 1.67                  | +8.00 to -17.00  | -6.00 | -17.00  | +0.50 to +1.25 | 70  | 3     |
| 1.74                  | +16.00 to -20.00 | -6.00 | -20.00  | +0.50 to +1.25 | 70  | 3     |
| PhotoFusion X         | SPH              | CYL   | COM RX. | ADD            | DIA | PRISM |
| 1.50                  | +5.00 to -7.00   | -4.00 | -7.00   | +0.50 to +1.25 | 75  | 3     |
| 1.60                  | +6.00 to -10.00  | -6.00 | -10.00  | +0.50 to +1.25 | 75  | 3     |
| 1.67                  | +8.00 to -12.00  | -6.00 | -12.00  | +0.50 to +1.25 | 75  | 3     |
| 1.74                  | +9.00 to -14.00  | -6.00 | -14.00  | +0.50 to +1.25 | 75  | 3     |
| AdaptiveSun           | SPH              | CYL   | COM RX. | ADD            | DIA | PRISM |
| 1.50                  | +5.00 to -7.00   | -4.00 | -7.00   | +0.50 to +1.25 | 75  | 3     |
| 1.60                  | +6.00 to -10.00  | -6.00 | -10.00  | +0.50 to +1.25 | 75  | 3     |
| 1.67                  | +8.00 to -12.00  | -6.00 | -12.00  | +0.50 to +1.25 | 70  | 3     |
| AdaptiveSun Polarized | SPH              | CYL   | COM RX. | ADD            | DIA | PRISM |
| 1.50                  | +5.00 to -8.00   | -4.00 | -8.00   | +0.50 to +1.25 | 80  | 3     |
| 1.60                  | +6.00 to -11.00  | -6.00 | -11.00  | +0.50 to +1.25 | 80  | 3     |
| Polarized             | SPH              | CYL   | COM RX. | ADD            | DIA | PRISM |
| 1.50                  | +5.00 to -8.00   | -4.00 | -8.00   | +0.50 to +1.25 | 80  | 3     |
| 1.60                  | +6.00 to -11.00  | -6.00 | -11.00  | +0.50 to +1.25 | 80  | 3     |
| 1.67                  | +8.00 to -12.00  | -6.00 | -12.00  | +0.50 to +1.25 | 80  | 3     |

\* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

\* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Price quoted is as per pair in Ringgit Malaysia and subject to prevailing tax. Price is subject to change without prior notice.

| ZEISS ClearMind<br>Progressive Individual 3 | BlueGuard |      |      |      |                    |      |      |        |
|---------------------------------------------|-----------|------|------|------|--------------------|------|------|--------|
|                                             | 1.50      | 1.60 | 1.67 | 1.74 |                    |      |      |        |
| DuraVision® Plus Gold UV                    | 5238      | 5788 | 6538 | 7238 |                    |      |      |        |
| DuraVision® Plus Platinum UV                | 5050      | 5600 | 6350 | 7050 | Now with BlueGuard |      |      |        |
| Clear   PhotoFusion X   Tinted              | UVProtect |      |      |      | PhotoFusion X      |      |      |        |
|                                             | 1.50      | 1.60 | 1.67 | 1.74 | 1.50               | 1.60 | 1.67 | 1.74** |
| DuraVision® Plus Mirror UV                  | 5150      | 5700 | -    | -    | -                  | -    | -    | -      |
| DuraVision® Plus Sun UV                     | 5050      | 5600 | 6350 | -    | -                  | -    | -    | -      |
| DuraVision® Plus Flash UV                   | 5150      | 5700 | 6450 | -    | 5650               | 6200 | 6950 | -      |
| DuraVision® Plus Gold UV                    | 5150      | 5700 | 6450 | 7150 | 5750               | 6300 | 7050 | 7750   |
| DuraVision® Plus Platinum UV                | 4950      | 5500 | 6250 | 6950 | 5550               | 6100 | 6850 | 7550   |

\* List price for DuraVision® Plus Flash, DuraVision® Plus Sun, DuraVision® Plus Mirror on UV Protect is inclusive of Standard Tint colour, excluding Syklet and Special Tint.

\* Top-up RM100/pair for Skylet (Fun, Road and Sport) & Special Tint (ProGolf, ProGolf Gradient and Double Gradient Tints)

\*\*PhotoFusion X 1.74 available in Grey, Extra Grey, Brown and Black only.

| AdaptiveSun                  | AdaptiveSun |      |      |      | AdaptiveSun Polarized |      |      |
|------------------------------|-------------|------|------|------|-----------------------|------|------|
|                              | 1.50        | 1.60 | 1.67 | 1.74 | 1.50                  | 1.60 | 1.67 |
| DuraVision® Plus Mirror UV   | 5850        | 6400 | -    | -    | 6000                  | 6550 | -    |
| DuraVision® Plus Sun UV      | 5750        | 6300 | 7050 | -    | 5900                  | 6450 | -    |
| DuraVision® Plus Flash UV    | 5850        | 6400 | 7150 | -    | 6000                  | 6550 | -    |
| DuraVision® Plus Gold UV     | 5950        | 6500 | 7250 | -    | 6100                  | 6650 | -    |
| DuraVision® Plus Platinum UV | 5750        | 6300 | 7050 | -    | 5900                  | 6450 | -    |
| Polarized                    | 1.50        | 1.60 | 1.67 | 1.74 |                       |      |      |
| DuraVision® Plus Mirror UV   | 5650        | 6200 | -    | -    |                       |      |      |
| DuraVision® Plus Sun UV      | 5550        | 6100 | 6850 | -    |                       |      |      |
| DuraVision® Plus Flash UV    | 5650        | 6200 | 6950 | -    |                       |      |      |
| DuraVision® Plus Gold UV     | 5750        | 6300 | 7050 | -    |                       |      |      |
| DuraVision® Plus Platinum UV | 5550        | 6100 | 6850 | -    |                       |      |      |

\* Require customer's date of birth (e.g. 14 Feb 1970), Frame Shape and Size, Pantoscopic angle, Vertex Distance, Bow Angle, FrameFit value (FFV from -1.0 to 6.0) for lens optimization.

\* Require customer's Vision Profile: Balance, Intermediate or Near for lifestyle optimization.

\* AdaptiveSun Blue is only available in 1.5 and 1.6 index.

\* AdaptiveSun and AdaptiveSun Polarized is not available in DuraVision® Plus Mirror UV Strong Blue.

\* Tinted colour is required for DuraVision® Plus Mirror UV, DuraVision® Plus Sun UV, DuraVision® Plus Flash UV (except for PhotoFusion X)

| BlueGuard / UVProtect | SPH              | CYL   | COM RX. | ADD            | DIA | PRISM |
|-----------------------|------------------|-------|---------|----------------|-----|-------|
| 1.50                  | +6.00 to -7.00   | -6.00 | -7.00   | +0.75 to +3.50 | 75  | 3     |
| 1.60                  | +10.00 to -10.00 | -6.00 | -10.00  | +0.75 to +4.00 | 75  | 3     |
| 1.67                  | +8.00 to -17.00  | -6.00 | -17.00  | +0.75 to +3.50 | 70  | 3     |
| 1.74                  | +16.00 to -20.00 | -6.00 | -20.00  | +0.75 to +3.50 | 70  | 3     |
| PhotoFusion X         | SPH              | CYL   | COM RX. | ADD            | DIA | PRISM |
| 1.50                  | +5.00 to -7.00   | -4.00 | -7.00   | +0.75 to +3.50 | 75  | 3     |
| 1.60                  | +6.00 to -10.00  | -6.00 | -10.00  | +0.75 to +4.00 | 75  | 3     |
| 1.67                  | +8.00 to -12.00  | -6.00 | -12.00  | +0.75 to +3.50 | 75  | 3     |
| 1.74                  | +9.00 to -14.00  | -6.00 | -14.00  | +0.75 to +3.50 | 75  | 3     |
| AdaptiveSun           | SPH              | CYL   | COM RX. | ADD            | DIA | PRISM |
| 1.50                  | +5.00 to -7.00   | -4.00 | -7.00   | +0.75 to +3.50 | 75  | 3     |
| 1.60                  | +6.00 to -10.00  | -6.00 | -10.00  | +0.75 to +4.00 | 75  | 3     |
| 1.67                  | +8.00 to -12.00  | -6.00 | -12.00  | +0.75 to +3.50 | 70  | 3     |
| AdaptiveSun Polarized | SPH              | CYL   | COM RX. | ADD            | DIA | PRISM |
| 1.50                  | +5.00 to -8.00   | -4.00 | -8.00   | +0.75 to +3.50 | 80  | 3     |
| 1.60                  | +6.00 to -11.00  | -6.00 | -11.00  | +0.75 to +3.50 | 80  | 3     |
| Polarized             | SPH              | CYL   | COM RX. | ADD            | DIA | PRISM |
| 1.50                  | +5.00 to -8.00   | -4.00 | -8.00   | +0.75 to +3.50 | 80  | 3     |
| 1.60                  | +6.00 to -11.00  | -6.00 | -11.00  | +0.75 to +3.50 | 80  | 3     |
| 1.67                  | +8.00 to -12.00  | -6.00 | -12.00  | +0.75 to +3.50 | 80  | 3     |

\* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

\* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Price quoted is as per pair in Ringgit Malaysia and subject to prevailing tax. Price is subject to change without prior notice.

| ZEISS ClearMind Progressive<br>Individual 3 Mineral Lens | Clear |      | Umbramatic |
|----------------------------------------------------------|-------|------|------------|
|                                                          | 1.60  | 1.80 | 1.60       |
| Super ET                                                 | 5600  | 7050 | 6200       |

\* Require customer's date of birth (e.g. 14 Feb 1970), Frame Shape and Size, Pantoscopic Angle, Back Vertex Distance, FrameFit (FFV: 0.0 - 6.0) for lens optimization.

\* Require 14 days delivery.

\* Tinted colour is unavailable.

| Clear      | SPH             | CYL   | COM RX. | ADD            | DIA | PRISM |
|------------|-----------------|-------|---------|----------------|-----|-------|
| 1.60       | +6.50 to -11.00 | -6.00 | -11.00  | +0.75 to +3.50 | 70  | 3     |
| 1.80       | 0.00 to -20.00  | -8.00 | -20.00  | +0.75 to +3.50 | 65  | 3     |
| Umbramatic | SPH             | CYL   | COM RX. | ADD            | DIA | PRISM |
| 1.60       | +6.50 to -9.00  | -4.00 | -9.00   | +0.75 to +3.50 | 70  | 3     |

\* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

You are advised to use VISUSTORE® to verify the availability of maximum diameter.

| ZEISS ClearMind Digital<br>Individual 3 Mineral Lens | Clear |      | Umbramatic |
|------------------------------------------------------|-------|------|------------|
|                                                      | 1.60  | 1.80 | 1.60       |
| Super ET                                             | 2190  | 3070 | 2790       |

\* Require customer's date of birth (e.g. 14 Feb 1970), Frame Shape and Size, Pantoscopic Angle, Back Vertex Distance, FrameFit (FFV: 0.0 - 6.0) for lens optimization.

\* Require 14 days delivery.


\* Tinted colour is unavailable.

| Clear      | SPH             | CYL   | COM RX. | ADD            | DIA | PRISM |
|------------|-----------------|-------|---------|----------------|-----|-------|
| 1.60       | +6.50 to -11.00 | -6.00 | -11.00  | +0.50 to +1.25 | 70  | 3     |
| 1.80       | 0.00 to -20.00  | -8.00 | -20.00  | +0.50 to +1.25 | 65  | 3     |
| Umbramatic | SPH             | CYL   | COM RX. | ADD            | DIA | PRISM |
| 1.60       | +6.50 to -9.00  | -4.00 | -9.00   | +0.50 to +1.25 | 70  | 3     |


\* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

You are advised to use VISUSTORE® to verify the availability of maximum diameter.


# Stamp and Engravings



ZEISS ClearMind Single Vision lenses




ZEISS Digital ClearMind Lenses lenses




ZEISS Digital ClearMind Individual 3 Lenses

# Stamp and Engravings



ZEISS Progressive ClearMind Individual 3 lens



ZEISS Progressive ClearMind Pure lens

ZEISS Progressive ClearMind Plus / Superb lens

## ZEISS ClearMind Lenses Lens Engraving

| Nasal Engraving                            | ClearMind |
|--------------------------------------------|-----------|
| ZEISS Single Vision ClearMind              | VCMxx     |
| ZEISS Single Vision ClearMind Individual 3 | VCIdx     |
| ZEISS Digital ClearMind                    | BCMxx     |
| ZEISS Digital Asiana* ClearMind            | DCAxx     |
| ZEISS Digital ClearMind Individual 3       | BCIdx     |
| ZEISS Progressive ClearMind Pure           | CPUxx     |
| ZEISS Progressive ClearMind Plus           | CMPxx     |
| ZEISS Progressive ClearMind Plus Asiana*   | CPAxx     |
| ZEISS Progressive ClearMind Superb         | CMSxx     |
| ZEISS Progressive ClearMind Individual 3   | CMDxx     |

\*Asiana version only available for some countries.