

Natural vision, fast adaptation.

ZEISS Light 2 Lenses

ZEISS

Seeing beyond

Light 2
Lenses



**Your light way to progressive lenses -
for the first timers and mature wearers.**

- Natural vision at all distances – near and far
- Fast adaptation
- Lightweight and thin lenses for comfort at all times

ZEISS Light 2 lenses

Quality optics with a focus on life

ZEISS Light 2 lenses are designed for your customers to switch easily to the progressive lens category while restoring their near and far vision effortlessly. It is ideal for first-timers or even mature wearers. It only takes one pair of glasses, and it won't cost a fortune. **Let your customers experience the world, entirely.**

Natural vision, easy adaptation

ZEISS Progressive Light 2 lenses have a moderate-soft design, which means they have a medium field of clear view but show a soft increase of blur distributed across the peripheral areas. As a result, the transition to blur is moderate-soft. **This allows for fast adaptation and natural vision at all distances for the wearer.**

Our promise

Our top priority and focus is to meet your customers' expectations and optimally respond to their ever-evolving needs. The engraving on each **ZEISS lens guarantees a high level of quality and symbolises** our legacy of precision in optics.

Great looks

ZEISS' commitment to optimizing lens aesthetics ensures that your customer's lenses are **lightweight and thin** for comfort at all times in a lens that looks good.

90%

of wearers had a pleasant **feel-good experience from the very first time** they put on their ZEISS Light 2 lenses.²



Target wearer

ZEISS Digital Light 2 lenses

For people who need easy reading of digital devices (smartphones, tablets etc.)

NEW



ZEISS Progressive Light 2 Lenses

For people experiencing presbyopia

UPGRADED




² Survey among purchasers of ZEISS Light 2 lenses in India (N=195, Digital and Progressive lenses), 2023, Carl Zeiss Vision International GmbH, DE (unpublished, data on file).

ZEISS EasyFocus technology

Innovative technologies for mature eyes

Our signature ZEISS EasyFocus technology considers the eye's inclination, the power profile and the binocular vision of the left and right eyes (of a progressive lens wearer) to optimise the lens for less blur in the periphery.

In addition, it is built on ZEISS' know-how in optics of the lens-eye-system, the age-related vision needs of first-time and mature progressive lens wearer and lens aesthetics optimisation.




The best of two worlds: medium-soft design fingerprint

ZEISS Light 2 lenses are a light way to progressive lenses because its medium-soft design fingerprint combines the benefits of hard and soft design of progressive lenses:

- ✓ Distance zones are wide and clear
- ✓ Very low peripheral distortion
- ✓ Enlarged near zones
- ✓ Accurately shaped to wearer's needs

up to **69%**
larger binocular vision zones than
our standard progressive lenses. ¹



1. Survey among eye care professionals, who sold ZEISS Light 2 lenses in India (N=42), 2023, Carl Zeiss Vision International GmbH, DE (unpublished, data on file).

ZEISS Progressive Light 2 Lenses



Seeing beyond



ZEISS Progressive Light 2 Lenses

Target group

ZEISS Progressive Light 2 lenses are designed to meet the vision needs of spectacle lens wearers who are experiencing presbyopia and need an addition power of more than 0.75D.

Benefits

- Restores vision effortlessly without costing a fortune.
- Simply one pair of glasses for everyday wear.
- **Fast adaptation.**
- **Natural and clear vision at all distances.**
- Easy reading of digital devices.
- **Lightweight and thin lenses** for comfort at all times.

Advantages

- ZEISS Progressive Light 2 lenses have a moderate-soft design, which means they have a medium field of clear view but show a soft increase in blur distributed across the peripheral areas. As a result, the transition to blur is moderate-soft. This allows for fast adaptation and clear vision at any distance for the wearer.
- ZEISS Progressive Light 2 lenses start from +0.75D up to +4.00D addition power and are available in three variations to meet the different visual needs and purchasing power of the wearer.



Large clear zones



Low blur in lens periphery



Accommodative support for comfortable reading with large near zone.

94%

of wearers found it **easy to switch** from their previous lens solution.²

68%

of wearers **adapted within** a few hours.²



2. Survey among purchasers of ZEISS Light 2 lenses in India (N=195, Digital and Progressive lenses), 2023, Carl Zeiss Vision International GmbH, DE (unpublished, data on file).

ZEISS Digital Light 2 Lenses



Seeing beyond



ZEISS Digital Light 2 Lenses

Natural vision, fast adaptation - first-timers.

Target group

ZEISS Digital Light 2 lenses are designed to meet the vision needs of spectacle lens wearers who experience tired of eyes and blurred near vision at the end of the day.

Benefits

- Restores vision effortlessly without costing a fortune.
- Simply one pair of glasses for everyday wear.
- Natural and clear vision at all distances.
- Easy reading of digital devices (smartphones, tablets etc.)
- Fast adaptation - it's even easy for single vision or non-lens wearers to transition into progressive wear.
- Lightweight and thin lenses for comfort at all times.

Advantages

- ZEISS Digital Light 2 lenses provide accommodative support in the lower lens part with a fixed addition of +0.50D or +1.00D for comfortable reading of close-up objects (specifically digital devices).
- The lens design is optimised for non-lens or single vision lens wearers whose eyes are unaccustomed to different visual zones. This allows for fast adaptation and clear vision at any distance for the wearer.



Large clear zones



Low blur in lens periphery



Accommodative support for comfortable reading with large near zone.

97%

of wearers confirm
**the lenses are
lightweight and
thing for comfort
at all times.¹**

96%

of wearers agree
**that their vision
feels natural.¹**

up
to **54%**

**larger binocular
vision zones** than
our standard
progressive lenses.²



1. Survey among purchasers of ZEISS Light 2 lenses in India (N=195, Digital & Progressive lenses), 2023, Carl Zeiss Vision International GmbH, DE (unpublished, data on file).

2. Calculation and Analysis by Technology and Innovation, Carl Zeiss Vision International GmbH, DE 2023. ZEISS Progressive and Digital Light 2 3Dve lenses compared to synchrony Ultra HD, synchrony Easy HD, ZEISS Progressive GT2 semi-finished and ZEISS Progressive Choice lenses, (unpublished, data on file).

Natural vision, fast adaptation.

The best ZEISS Light 2
lenses for you.



Natural and clear vision
at all distances



Fast adaptation



Lightweight and thin lenses
for comfort at all times



92% of wearers had a pleasant **feel-good experience**
from the first time they put on the lenses.¹

You are	For people who need easy reading of digital devices (smartphones, tablets etc.)		Designed for progressive lens wearers with $\geq +0.75\text{D}$ addition power (up to $+4.00\text{D}$ addition power)	
	OUR GOOD SOLUTION.	OUR BETTER SOLUTION.	OUR GOOD SOLUTION.	OUR BASIC SOLUTION.
Lens options	ZEISS Digital Light 2 3D lenses	ZEISS Progressive Light 2 3Dv lenses	ZEISS Progressive Light 2 3D lenses	ZEISS Progressive Light 2 D lenses
Vision zones	- Large clear zones - Accommodative support in the lower lens part for relaxed reading - Soft peripheral blur zones - Optimised for digital device use	- Larger clear zones than ZEISS Progressive Light 2 3D when fitted in near and intermediate zone - Soft peripheral blur - Optimised for digital device use	- Larger clear zones than ZEISS Progressive Light 2 D when fitted in near vision zone - Soft peripheral blur - Optimised for digital device use	- Medium clear zones - Soft peripheral blur
Frame choice	 For standard frames, available in 3 fitting heights (S, M and L).	 For any frame of your choice, fitted in near and intermediate zones.	 For standard frames, available in 3 fitting heights (S, M and L).	 For standard frames, available in 2 fitting heights (S and L).

1. Survey among purchasers of ZEISS Light 2 lenses in India (N=195, Digital & Progressive lenses), 2023, Carl Zeiss Vision International GmbH, DE (unpublished, data on file)

ZEISS Progressive Light 2 3Dv	BlueGuard							
	1.50	1.60	1.67	1.74				
DuraVision® Plus Gold UV	1900	2300	2660	3050				
DuraVision® Plus Platinum UV	1700	2100	2460	2850				
DuraVision® Plus Chrome UV	1600	2000	2360	2750				Now with BlueGuard
Clear PhotoFusion X Tinted	UVProtect				PhotoFusion X			
	1.50	1.60	1.67	1.74	1.50	1.60	1.67	1.74**
DuraVision® Plus Mirror UV	1800	2200	-	-	-	-	-	-
DuraVision® Plus Sun UV	1700	2100	2460	-	-	-	-	-
DuraVision® Plus Flash UV	1800	2200	2560	-	2300	2700	3060	-
DuraVision® Plus Gold UV	1800	2200	2560	2950	2400	2800	3160	3550
DuraVision® Plus Platinum UV	1600	2000	2360	2750	2200	2600	2960	3350
DuraVision® Plus Chrome UV	1500	1900	2260	2650	2100	2500	2860	3250

* List price for DuraVision® Plus Flash, DuraVision® Plus Sun, DuraVision® Plus Mirror on UV Protect is inclusive of Standard Tint colour, excluding Syklet and Special Tint.

* Top-up RM100/pair for Skylet (Fun, Road and Sport) & Special Tint (ProGolf, ProGolf Gradient and Double Gradient Tints)

**PhotoFusion X 1.74 available in Grey, Extra Grey, Brown and Black only.

AdaptiveSun	AdaptiveSun				AdaptiveSun Polarized		
	1.50	1.60	1.67	1.74	1.50	1.60	1.67
DuraVision® Plus Mirror UV	2500	2900	-	-	2650	3050	-
DuraVision® Plus Sun UV	2400	2800	3160	-	2550	2950	-
DuraVision® Plus Flash UV	2500	2900	3260	-	2650	3050	-
DuraVision® Plus Gold UV	2600	3000	3360	-	2750	3150	-
DuraVision® Plus Platinum UV	2400	2800	3160	-	2550	2950	-
DuraVision® Plus Chrome UV	2300	2700	3060	-	2450	2850	
Polarized	1.50	1.60	1.67	1.74			
DuraVision® Plus Mirror UV	2300	2700	-	-			
DuraVision® Plus Sun UV	2200	2600	2960	-			
DuraVision® Plus Flash UV	2300	2700	3060	-			
DuraVision® Plus Gold UV	2400	2800	3160	-			
DuraVision® Plus Platinum UV	2200	2600	2960	-			
DuraVision® Plus Chrome UV	2100	2500	2860	-			

* Require customer's date of birth (e.g. 14 Feb 1996), Frame Shape and Size, for lens optimization

* Require Frame Shape and Size and FrameFit value (FFV from -1.0 to 6.0) for lens optimization

* Tinted colour is required for DuraVision® Plus Mirror UV, DuraVision® Plus Sun UV, DuraVision® Plus Flash UV (except for PhotoFusion X)

BlueGuard / UVProtect	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	70	3
1.74	+9.00 to -14.00	-6.00	-14.00	+0.75 to +3.50	75	3
PhotoFusion X	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	75	3
1.74	+9.00 to -14.00	-6.00	-14.00	+0.75 to +3.50	75	3
AdaptiveSun	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	70	3
AdaptiveSun Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.75 to +3.50	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.75 to +3.50	80	3
Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.75 to +3.50	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.75 to +3.50	80	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	80	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Price quoted is as per pair in Ringgit Malaysia and subject to prevailing tax. Price is subject to change without prior notice.

ZEISS Progressive Light 2 3D	BlueGuard							
	1.50	1.60	1.67	1.74				
DuraVision® Plus Gold UV	1420	1820	2160	2750				
DuraVision® Plus Platinum UV	1220	1620	1960	2550				
DuraVision® Plus Chrome UV	1120	1520	1860	2450			Now with BlueGuard	
Clear PhotoFusion X Tinted	UVProtect				PhotoFusion X			
	1.50	1.60	1.67	1.74	1.50	1.60	1.67	1.74**
DuraVision® Plus Mirror UV	1320	1720	-	-	-	-	-	-
DuraVision® Plus Sun UV	1220	1620	1960	-	-	-	-	-
DuraVision® Plus Flash UV	1320	1720	2060	-	1820	2220	2560	-
DuraVision® Plus Gold UV	1320	1720	2060	2650	1920	2320	2660	3250
DuraVision® Plus Platinum UV	1120	1520	1860	2450	1720	2120	2460	3050
DuraVision® Plus Chrome UV	1020	1420	1760	2350	1620	2020	2360	2950

* List price for DuraVision® Plus Flash, DuraVision® Plus Sun, DuraVision® Plus Mirror on UV Protect is inclusive of Standard Tint colour, excluding Syklet and Special Tint.

* Top-up RM100/pair for Skylet (Fun, Road and Sport) & Special Tint (ProGolf, ProGolf Gradient and Double Gradient Tints)

**PhotoFusion X 1.74 available in Grey, Extra Grey, Brown and Black only.

AdaptiveSun	AdaptiveSun				AdaptiveSun Polarized		
	1.50	1.60	1.67	1.74	1.50	1.60	1.67
DuraVision® Plus Mirror UV	2020	2420	-	-	2170	2570	-
DuraVision® Plus Sun UV	1920	2320	2660	-	2070	2470	-
DuraVision® Plus Flash UV	2020	2420	2760	-	2170	2570	-
DuraVision® Plus Gold UV	2120	2520	2860	-	2270	2670	-
DuraVision® Plus Platinum UV	1920	2320	2660	-	2070	2470	-
DuraVision® Plus Chrome UV	1820	2220	2560	-	1970	2370	
Polarized	1.50	1.60	1.67	1.74			
DuraVision® Plus Mirror UV	1820	2220	-	-			
DuraVision® Plus Sun UV	1720	2120	2460	-			
DuraVision® Plus Flash UV	1820	2220	2560	-			
DuraVision® Plus Gold UV	1920	2320	2660	-			
DuraVision® Plus Platinum UV	1720	2120	2460	-			
DuraVision® Plus Chrome UV	1620	2020	2360	-			

* Require customer's date of birth (e.g. 14 Feb 1996), Frame Shape and Size for lens optimization

* Require specific lens corridor length: Long-14mm (min.fitting 18mm), Medium-12mm (min.fitting 16mm) or Short-10mm (min.fitting 14mm)

* Tinted colour is required for DuraVision® Plus Mirror UV, DuraVision® Plus Sun UV, DuraVision® Plus Flash UV (except for PhotoFusion X)

BlueGuard / UVProtect	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	70	3
1.74	+9.00 to -14.00	-6.00	-14.00	+0.75 to +3.50	75	3
PhotoFusion X	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	75	3
1.74	+9.00 to -14.00	-6.00	-14.00	+0.75 to +3.50	75	3
AdaptiveSun	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	70	3
AdaptiveSun Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.75 to +3.50	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.75 to +3.50	80	3
Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.75 to +3.50	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.75 to +3.50	80	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	80	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Price quoted is as per pair in Ringgit Malaysia and subject to prevailing tax. Price is subject to change without prior notice.

ZEISS Progressive Light 2 D	BlueGuard			Now with BlueGuard			
	1.50	1.60	1.67				
DuraVision® Plus Gold UV	1110	1490	1860				
DuraVision® Plus Platinum UV	910	1290	1660				
DuraVision® Plus Chrome UV	810	1190	1560				
Clear PhotoFusion X Tinted	UVProtect				PhotoFusion X		
	1.50	1.60	1.67		1.50	1.60	1.67
DuraVision® Plus Mirror UV	1010	1390	-		-	-	-
DuraVision® Plus Sun UV	910	1290	1660		-	-	-
DuraVision® Plus Flash UV	1010	1390	1760		1510	1890	2260
DuraVision® Plus Gold UV	1010	1390	1760	1610	1990	2360	
DuraVision® Plus Platinum UV	810	1190	1560	1410	1790	2160	
DuraVision® Plus Chrome UV	710	1090	1460	1310	1690	2060	

* Require customer's date of birth (e.g. 14 Feb 1996), Frame Shape and Size for lens optimization

* Require specific lens corridor length: Long-14mm (min.fitting 18mm) or Short-10mm (min.fitting 14mm)

* Tinted colour is required for DuraVision® Plus Mirror UV, DuraVision® Plus Sun UV, DuraVision® Plus Flash UV (except for PhotoFusion X)

* List price for DuraVision® Plus Flash, DuraVision® Plus Sun, DuraVision® Plus Mirror on UV Protect is inclusive of Standard Tint colour, excluding Syklet and Special Tint

* Top-up RM100/pair for Skylet (Skylet fun, SkyPol fun, Skylet Sport, SkyPol Sport) & Special Tint (ProGolf, ProGolf Gradient, Double Gradient Tints)

BlueGuard / UVProtect	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	70	3
PhotoFusion X	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	75	3
1.74	+9.00 to -14.00	-6.00	-14.00	+0.75 to +3.50	75	3
Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.75 to +3.50	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.75 to +3.50	80	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	80	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

ZEISS Digital Light 2 3D	BlueGuard							
	1.50	1.60	1.67	1.74				
DuraVision® Plus Gold UV	1000	1400	1600	2250				
DuraVision® Plus Platinum UV	800	1200	1400	2050				
DuraVision® Plus Chrome UV	700	1100	1300	1950				
					Now with BlueGuard			
Clear PhotoFusion X Tinted	UVProtect				PhotoFusion X			
	1.50	1.60	1.67	1.74	1.50	1.60	1.67	1.74**
DuraVision® Plus Mirror UV	900	1300	-	-	-	-	-	-
DuraVision® Plus Sun UV	800	1200	1400	-	-	-	-	-
DuraVision® Plus Flash UV	900	1300	1500	-	1400	1800	2000	-
DuraVision® Plus Gold UV	900	1300	1500	2150	1500	1900	2100	2750
DuraVision® Plus Platinum UV	700	1100	1300	1950	1300	1700	1900	2550
DuraVision® Plus Chrome UV	600	1000	1200	1850	1200	1600	1800	2450

* List price for DuraVision® Plus Flash, DuraVision® Plus Sun, DuraVision® Plus Mirror on UV Protect is inclusive of Standard Tint colour, excluding Syklet and Special Tint.

** Top-up RM100/pair for Skylet (Fun, Road and Sport) & Special Tint (ProGolf, ProGolf Gradient and Double Gradient Tints)

**PhotoFusion X 1.74 available in Grey, Extra Grey, Brown and Black only.

AdaptiveSun	AdaptiveSun				AdaptiveSun Polarized		
	1.50	1.60	1.67	1.74	1.50	1.60	1.67
DuraVision® Plus Mirror UV	1600	2000	-	-	1750	2150	-
DuraVision® Plus Sun UV	1500	1900	2100	-	1650	2050	-
DuraVision® Plus Flash UV	1600	2000	2200	-	1750	2150	-
DuraVision® Plus Gold UV	1700	2100	2300	-	1850	2250	-
DuraVision® Plus Platinum UV	1500	1900	2100	-	1650	2050	-
DuraVision® Plus Chrome UV	1400	1800	2000	-	1550	1950	-
Polarized	1.50	1.60	1.67	1.74			
DuraVision® Plus Mirror UV	1400	1800	-	-			
DuraVision® Plus Sun UV	1300	1700	1900	-			
DuraVision® Plus Flash UV	1400	1800	2000	-			
DuraVision® Plus Gold UV	1500	1900	2100	-			
DuraVision® Plus Platinum UV	1300	1700	1900	-			
DuraVision® Plus Chrome UV	1200	1600	1800	-			

* Require customer's date of birth (e.g. 14 Feb 1996), Frame Shape and Size for lens optimization

* Require specific lens corridor length: Long-14mm (min.fitting 18mm), Medium-12mm (min.fitting 16mm) or Short-10mm (min.fitting 14mm)

* Tinted colour is required for DuraVision® Plus Mirror UV, DuraVision® Plus Sun UV, DuraVision® Plus Flash UV (except for PhotoFusion X)

BlueGuard / UVProtect	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.50 / +1.00	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.50 / +1.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.50 / +1.00	70	3
1.74	+9.00 to -14.00	-6.00	-14.00	+0.50 / +1.00	70	3
PhotoFusion X	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.50 / +1.00	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.50 / +1.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.50 / +1.00	75	3
1.74	+9.00 to -14.00	-6.00	-14.00	+0.50 / +1.00	75	3
AdaptiveSun	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.50 / +1.00	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.50 / +1.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.50 / +1.00	70	3
AdaptiveSun Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.50 / +1.00	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.50 / +1.00	80	3
Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.50 / +1.00	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.50 / +1.00	80	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.50 / +1.00	80	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Price quoted is as per pair in Ringgit Malaysia and subject to prevailing tax. Price is subject to change without prior notice.