

ZEISS SmartLife PRO


Progressive Lenses



ZEISS SmartLife PRO Progressive Lens Portfolio

Scientifically proven to enable peripheral vision in natural dynamic interaction.

ZEISS SmartLife PRO Progressive Lenses are designed to meet the vision needs of people with presbyopia with a connected and on-the-move lifestyle. The optical performance in the lens periphery is designed for frequent changes of head and eye position driven by how people interact with their handheld devices. This new design fingerprint provides a smoother transition into the lens periphery with less perceived blur to enable peripheral vision in natural dynamic interaction.




Smoother transitions

Smoother transition into areas with more blur compared to current ZEISS Precision Progressive lenses.


Lower blur level

Decreased blur level for intermediate to near distances compared to ZEISS Precision lenses.

Standard Progressive Lens



ZEISS SmartLife PRO Progressive Lens



Consumer benefits


- High performance in dynamic vision in the periphery.
- 9 out of 10 rated the quality of vision with ZEISS SmartLife lenses positive.
- 4 out of 5 experienced smooth vision from near to far across all viewing zones.
- 8 of out 10 consumers adapted very fast to their new lenses, (within 1 day).



ZEISS SmartLife PRO Progressive Pure



**ZEISS SmartLife PRO
Progressive Pure**
Optimised for the eyes



Digital Inside® Technology
Pure vision with enhanced performance
for the digital world.



One of the main trends in today's world is the increased use of digital devices. Interaction with computers, mobile phones and tablets has become embedded in our business and personal lives over the past few years.

Human eyes are challenged like never before

The reading distance when using digital devices is closer than with books or any other print media. Eye movements are fast and dynamic, from near to far and back, all day long. This must be taken into consideration in order to avoid severe eye strain and fatigue.

What is the problem with conventional progressive lenses in today's world?



Progressive lens optimised for conventional reading behaviour:

Near vision zone and inset adjusted for reading distance with print media.



Progressive lens not optimised for digital reading behaviour:

Near vision zone is not compatible with the shorter reading distance needed for digital devices. The result: an unnatural head and body posture is required to move digital devices into the sharp and clear near vision zone.

ZEISS SmartLife Progressive Pure responds to this new challenge with enhanced vision for the digital world



Digital Inside™ Technology Progressive lens optimised for conventional and digital reading behaviour:

Near vision zone is extended vertically and horizontally for both reading distances – conventional print media and digital devices. The result : relaxed and natural head and body posture.