

Physiologically fitted to the wearer's face.

ZEISS Individual Lenses



Individual
Lenses

Consumer benefits



Optimised binocular vision zones and full utilisation of 3D vision potential




Excellent vision performance for using digital devices: less eye strain, comfortable vision and easy focusing on near zones



Wide selection of frames

Progressive SmartLife Individual

Optimised for the eyes + frame + face + main daily activities



IndividualFit™ & Luminance Design® Technology

Tailor-made to fit all individual requirements.



Among the many influencing trends today, more and more people are enjoying an individualised way of life. And they are asking for the appropriate tailor-made solutions.

IndividualFit™ Technology

IndividualFit™ Technology takes into account the individual's main daily activities and provides a choice of three lens designs to suit the wearer's needs.



Near:

Optimised design for prolonged near vision activities.

- Optimum near vision zone
- Wide distance vision zone maintained
- Good intermediate vision zone



Intermediate:

Optimised design for dynamic and intermediate activities.

- Optimum intermediate vision zone
- Wide distance vision zone maintained
- Good near vision zone



Balanced:

Optimised design for all-round activities.

- Balanced vision ranges
- Wide distance, intermediate and near vision zones

In addition, ZEISS Progressive SmartLife Individual incorporates the individually measured reading distance and the individual wrap angle within the lens design.

Luminance Design™ Technology


Knowing that pupil size changes in different light conditions, the new Luminance Design™ Technology factors in an average pupil size of bundle for more natural vision. 3.3 mm. It calculates the lens design bundle by



Small pupil during the day




Mid-size pupil in mesopic conditions



Large pupil at night



Seeing beyond



ZEISS SmartLife Individual 3

with Intelligence Augmented Design[®]
technology by ZEISS

ZEISS SmartLife Individual 3 lenses feature ZEISS Intelligence Augmented Design™ technology.

ZEISS SmartLife Individual 3 is the next generation of tailor-made lenses. The evolutionary design makes use of intelligence augmentation to take unique visual requirements, behavior and parameters into account in the highest level of detail. These lenses are so precisely customized – **they are simply the best ZEISS lenses ever.**

Our calculations and studies show:

- 85% of wearer confirmed that they have full freedom of vision in any distance & direction9.
- Up to 29% better optical performance than previous SmartLife Individual lenses.
- 8 out of 10 wearers confirmed they had comfortable vision within minutes of wearing9.
- Up to 79% more suited to the visual needs and behavior of wearers compared to previous SmartLife Individual lenses.

ZEISS SmartLife Individual 3 lenses consider oculomotoric behavior together with the wearer’s individual parameters including anatomy, physiology as well as habituation, resulting in an ultra-refined lens.


The lenses make use of new **ZEISS Intelligence Augmented Design technology** to match an individual wearer’s data to a set of **12.5million data points of visual behavior** in specific age groups. This allow ZEISS to predict how an individual will use their eyes in their daily life and to **create their best matched lens design.**

Product availability

ZEISS SmartLife Individual 3 lenses are available for all wearers in SmartLife Single Vision, SmartLife Digital and SmartLife Progressive lenses regardless of material, and can be combined with all ZEISS DuraVision coatings, tints, and ZEISS PhotoFusion X.

Today


One Design - All Ages




Best balanced zone sizes for all age groups.

New


Best Match Design Age 47



Best Match Design Age 55



Best Match Design Age 73



Personalized zone sizes build on age-related differences in today's visual behavior.

ZEISS SmartLife Individual 3 lenses

The magic of individual lenses truly lies in the combination of technologies. We utilize up to 8 different technologies to create one highly refined ZEISS SmartLife Individual 3 lens. The complex interaction of the technology and the optimization process should not be seen in isolation - but rather as essential ingredients that result in the best possible product for the wearer.

FOCAL TYPE:		ZEISS Single Vision SmartLife Individual 3	ZEISS Digital SmartLife Individual 3	ZEISS Progressive SmartLife Individual 3
Technologies				
 ZEISS SmartView 2.0 technology Optimised for all-day visual comfort for a connected and on-the-move lifestyle. 				
 ZEISS Luminance Design 2.0 technology Takes average light conditions and age-related pupil diameter into account.				
 ZEISS Digital Inside technology Optimised near zone for better reading on digital devices.				
 ZEISS FrameFit+ technology Optimisation of the corridor for free frame choice.			 <i>variable corridor based on fitting height</i>	
 ZEISS Adaptation Control technology Optimised for faster adaptation from old to new lenses.				
 ZEISS FaceFit technology Fitted to the wearer's unique facial anatomy.				
 ZEISS IndividualFit technology Optimised for specific daily activities.				
 ZEISS Intelligence Augmented Design (IAD) technology  Optimised to wearer's personal visual behaviour and requirements.				

NEW

The paradigm shift of individual lenses

The paradigm shift of individual lenses

Brute force measurement makes way for data science

1


Measurement and individual data determined by ECP



Measurement and determination of individual data. Choose from 5 ZEISS SmartLife Individual Lenses (1xSV/ 1xDigi/ 3xPAL [B/N/I])

2


Individual data + ZEISS Intelligence Augmented Design technology + ZEISS FrameFit technology leads to target design.



The chosen SmartLife lens will be enriched with intelligence augmentation to generate a more precise design target.
→ "Individual" visual demand chosen from **1 out of 44.000** possible target designs

3

Individual data is used for lens performance optimisation



The lens performance will be optimised based on anatomical and frame data for wide fields of view
→ Individual lens performance based on anatomy and frame



12.5 mio
data points

of real-life habitual, lifestyle and oculomotoric behaviour were analysed and combined with anatomical and physiological information to create the new

ZEISS Intelligence Augmented Design® technology




ZEISS SmartLife Individual 3 Lenses

Technologies overview

Process step	Objectives	Supporting technologies
1 Determination of wearer's personal vision profile. 	The eyecare professional captures and sends all individual parameters to ZEISS together with the choice of the ZEISS SmartLife Individual lens. This can be a single vision, digital or a progressive individual lens.	ZEISS Adaptation Control technology (intergrated in the FrameFit calculator) will support the determination of the corridor length for our progressive & digital lenses.
2 Build on age-related differences in today's visual behaviour 	12.5mio data points of the real-life habitual, lifestyle & oculomotoric behaviour data were analysed and translated into age-dependent profiles as starting point for any ZEISS SmartLife Individual lens.	ZEISS Digital Inside technology & ZEISS SmartView 2.0 technology form the base of the design. It considers key elements of the visual behaviour and visual needs of people of any age in today's connected and on-the-move lifestyle and ZEISS IndividualFit technology is recommended for individual zone adaptation for wearers with specific visual tasks.
3 Create best match lens design target 	The age-dependent visual profiles will be combined with personal, anatomical and physiological information of the wearer to create the best matched lens design target out of 44k possibilities.	ZEISS Intelligence Augmented Design technology identifies the age-dependent visual behaviour profiles, check this against the optometric status and habituation of the individual wearer and employs ZEISS FrameFit+ technology to customise the lens corridor to the chosen frame for the best match design target.
4 Optimise final lens performance 	All the provided individual anatomical & frame data will be used to optimise the lens performance to fit the chosen target design at its best.	ZEISS Luminance Design 2.0 technology & ZEISS FaceFit technology optimise the real lens in an interactive process to come as close as possible to the chosen design target.
5 ZEISS SmartLife Individual 3 lens 	ZEISS SmartLife Individual 3 lens is as close as possible to the target design of the wearer to provide full freedom of clear vision in any distance & direction.	

Comfortable vision for busy eyes.

The best ZEISS SmartLife lenses for you.

	Lens options	Vision zones	Frame Choice
Designed for people with blurry near and far vision, and >1.25 ADD power. 45+	OUR BEST SOLUTION. ZEISS Progressive SmartLife Individual 3 lenses		 For any frame of your choice, adapted to your individual anatomy and special vision needs .
		■ Larger clear zones than ZEISS Progressive SmartLife Superb when fitted to individual anatomy & vision needs. ■ Smooth peripheral blur zones. ■ Optimised for digital device use.	
	OUR BETTER SOLUTION. ZEISS Progressive SmartLife Superb lenses		 For any frame of your choice, adapted to your individual anatomy .
		■ Larger clear zones than ZEISS Progressive SmartLife Plus when fitted to individual anatomy. ■ Smooth peripheral blur. ■ Optimised for digital device use.	
	OUR GOOD SOLUTION. ZEISS Progressive SmartLife Plus lenses		 For any frame of your choice, adapted to your individual anatomy .
		■ Larger clear zones than ZEISS Progressive SmartLife Pure when fitted in near vision zone. ■ Smooth peripheral blur. ■ Optimised for digital device use.	
	OUR BASIC SOLUTION. ZEISS Progressive SmartLife Pure lenses		 For standard frame available in 3 fitting heights (14/16/18 mm) .
		■ Large clear zones. ■ Smooth peripheral blur. ■ Optimised for digital device usage.	

1 Compared to ZEISS Progressive SmartLife Superb, (data on file).
The frames left to right: ZEISS Fall Winter 2023 collection: ZS23533; ZEISS Fall Winter 2023 collection: ZS23128; ZEISS Fall Winter 2022 collection: ZS22524; generic frame.

ZEISS Single Vision SmartLife Individual 3	BlueGuard							
	1.50	1.60	1.67	1.74				
DuraVision® Plus Gold UV	1780	2050	2270	2950				
DuraVision® Plus Platinum UV	1580	1850	2070	2750	Now with BlueGuard			
Clear PhotoFusion X Tinted	UVProtect				PhotoFusion X			
	1.50	1.60	1.67	1.74	1.50	1.60	1.67	1.74**
DuraVision® Plus Mirror UV	1680	1950	-	-	-	-	-	-
DuraVision® Plus Sun UV	1580	1850	2070	-	-	-	-	-
DuraVision® Plus Flash UV	1680	1950	2170	-	2180	2450	2670	-
DuraVision® Plus Gold UV	1680	1950	2170	2850	2280	2550	2770	3450
DuraVision® Plus Platinum UV	1480	1750	1970	2650	2080	2350	2570	3250

* List price for DuraVision Flash, DuraVision Sun, DuraVision Mirror on UV Protect is inclusive of Standard Tint colour, excluding Syklet and Special Tint.

* Top-up RM100/pair for Skylet (Fun, Road and Sport) & Special Tint (ProGolf, ProGolf Gradient and Double Gradient Tints)

**PhotoFusion X 1.74 available in Grey, Extra Grey, Brown and Black only.

AdaptiveSun	AdaptiveSun				AdaptiveSun Polarized		
	1.50	1.60	1.67	1.74	1.50	1.60	1.67
DuraVision® Plus Mirror UV	2380	2650	-	-	2530	2800	-
DuraVision® Plus Sun UV	2280	2550	2770	-	2430	2700	-
DuraVision® Plus Flash UV	2380	2650	2870	-	2530	2800	-
DuraVision® Plus Gold UV	2480	2750	2970	-	2630	2900	-
DuraVision® Plus Platinum UV	2280	2550	2770	-	2430	2700	-
Polarized	1.50	1.60	1.67	1.74			
DuraVision® Plus Mirror UV	2180	2450	-	-			
DuraVision® Plus Sun UV	2080	2350	2570	-			
DuraVision® Plus Flash UV	2180	2450	2670	-			
DuraVision® Plus Gold UV	2280	2550	2770	-			
DuraVision® Plus Platinum UV	2080	2350	2570	-			

* Require customer's date of birth (e.g. 14 Feb 1996), Frame Shape and Size, Pantoscopic angle, Wrap angle and Back Vertex Distance, for lens optimization

* AdaptiveSun Blue is only available in 1.5 and 1.6 index.

* AdaptiveSun and AdaptiveSun Polarized is not available in DuraVision Mirror UV Strong Blue.

* Tinted colour is required for DuraVision Mirror UV, DuraVision Sun UV, DuraVision Flash UV (except for PhotoFusion X)

BlueGuard / UVProtect	SPH	CYL	COM RX.	DIA	PRISM
1.50	+6.00 to -6.00	-4.00	-6.00	75	3
1.60	+10.00 to -10.00	-6.00	-10.00	75	3
1.67	+10.00 to -17.00	-6.00	-17.00	70	3
1.74	+16.00 to -20.00	-6.00	-20.00	70	3
PhotoFusion X	SPH	CYL	COM RX.	DIA	PRISM
1.50	+6.00 to -6.00	-4.00	-6.00	75	3
1.60	+8.00 to -10.00	-6.00	-10.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	75	3
1.74	+9.00 to -14.00	-6.00	-14.00	75	3
AdaptiveSun	SPH	CYL	COM RX.	DIA	PRISM
1.50	+6.00 to -6.00	-4.00	-6.00	75	3
1.60	+8.00 to -10.00	-6.00	-10.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	70	3
AdaptiveSun Polarized	SPH	CYL	COM RX.	DIA	PRISM
1.50	+6.00 to -6.00	-4.00	-6.00	80	3
1.60	+10.00 to -11.00	-6.00	-11.00	80	3
Polarized	SPH	CYL	COM RX.	DIA	PRISM
1.50	+8.00 to -8.00	-4.00	-8.00	80	3
1.60	+10.00 to -11.00	-6.00	-11.00	80	3
1.67	+11.00 to -12.00	-6.00	-12.00	80	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

ZEISS Digital SmartLife Individual 3	BlueGuard				Now with BlueGuard			
	1.50	1.60	1.67	1.74				
DuraVision® Plus Gold UV	1900	2200	2570	3000				
DuraVision® Plus Platinum UV	1700	2000	2370	2800				
Clear PhotoFusion X Tinted	UVProtect				PhotoFusion X			
	1.50	1.60	1.67	1.74	1.50	1.60	1.67	1.74**
DuraVision® Plus Mirror UV	1800	2100	-	-	-	-	-	-
DuraVision® Plus Sun UV	1700	2000	2370	-	-	-	-	-
DuraVision® Plus Flash UV	1800	2100	2470	-	2300	2600	2970	-
DuraVision® Plus Gold UV	1800	2100	2470	2900	2400	2700	3070	3500
DuraVision® Plus Platinum UV	1600	1900	2270	2700	2200	2500	2870	3300

* List price for DuraVision Flash, DuraVision Sun, DuraVision Mirror on UV Protect is inclusive of Standard Tint colour, excluding Syklet and Special Tint.

* Top-up RM100/pair for Skylet (Fun, Road and Sport) & Special Tint (ProGolf, ProGolf Gradient and Double Gradient Tints)

** PhotoFusion X 1.74 available in Grey, Extra Grey, Brown and Black only.

AdaptiveSun	AdaptiveSun				AdaptiveSun Polarized		
	1.50	1.60	1.67	1.74	1.50	1.60	1.67
DuraVision® Plus Mirror UV	2500	2800	-	-	2650	2950	-
DuraVision® Plus Sun UV	2400	2700	3070	-	2550	2850	-
DuraVision® Plus Flash UV	2500	2800	3170	-	2650	2950	-
DuraVision® Plus Gold UV	2600	2900	3270	-	2750	3050	-
DuraVision® Plus Platinum UV	2400	2700	3070	-	2550	2850	-
Polarized	1.50	1.60	1.67	1.74			
DuraVision® Plus Mirror UV	2300	2600	-	-			
DuraVision® Plus Sun UV	2200	2500	2870	-			
DuraVision® Plus Flash UV	2300	2600	2970	-			
DuraVision® Plus Gold UV	2400	2700	3070	-			
DuraVision® Plus Platinum UV	2200	2500	2870	-			

* Require customer's date of birth (e.g. 14 Feb 1996), Frame Shape and Size, Pantoscopic angle for lens optimization, FFV

* AdaptiveSun Blue is only available in 1.5 and 1.6 index.

* AdaptiveSun and AdaptiveSun Polarized is not available in DuraVision Mirror UV Strong Blue.

* Tinted colour is required for DuraVision Mirror UV, DuraVision Sun UV, DuraVision Flash UV (except for PhotoFusion X)

BlueGuard / UVProtect	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+6.00 to -7.00	-6.00	-7.00	+0.50 to +1.25	75	3
1.60	+10.00 to -10.00	-6.00	-10.00	+0.50 to +1.25	75	3
1.67	+8.00 to -17.00	-6.00	-17.00	+0.50 to +1.25	70	3
1.74	+16.00 to -20.00	-6.00	-20.00	+0.50 to +1.25	70	3
PhotoFusion X	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.50 to +1.25	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.50 to +1.25	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.50 to +1.25	75	3
1.74	+9.00 to -14.00	-6.00	-14.00	+0.50 to +1.25	75	3
AdaptiveSun	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.50 to +1.25	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.50 to +1.25	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.50 to +1.25	70	3
AdaptiveSun Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.50 to +1.25	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.50 to +1.25	80	3
Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.50 to +1.25	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.50 to +1.25	80	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.50 to +1.25	80	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Price quoted is as per pair in Ringgit Malaysia and subject to prevailing tax. Price is subject to change without prior notice.

ZEISS Progressive SmartLife Individual 3	BlueGuard							
	1.50	1.60	1.67	1.74				
DuraVision® Plus Gold UV	4850	5300	6000	6600				
DuraVision® Plus Platinum UV	4650	5100	5800	6400	Now with BlueGuard			
Clear PhotoFusion X Tinted	UVProtect				PhotoFusion X			
	1.50	1.60	1.67	1.74	1.50	1.60	1.67	1.74**
DuraVision® Plus Mirror UV	4750	5200	-	-	-	-	-	-
DuraVision® Plus Sun UV	4650	5100	5800	-	-	-	-	-
DuraVision® Plus Flash UV	4750	5200	5900	-	5250	5700	6400	-
DuraVision® Plus Gold UV	4750	5200	5900	6500	5350	5800	6500	7100
DuraVision® Plus Platinum UV	4550	5000	5700	6300	5150	5600	6300	6900

* List price for DuraVision Flash, DuraVision Sun, DuraVision Mirror on UV Protect is inclusive of Standard Tint colour, excluding Syklet and Special Tint.

* Top-up RM100/pair for Skylet (Fun, Road and Sport) & Special Tint (ProGolf, ProGolf Gradient and Double Gradient Tints)

** PhotoFusion X 1.74 available in Grey, Extra Grey, Brown and Black only.

AdaptiveSun	AdaptiveSun				AdaptiveSun Polarized		
	1.50	1.60	1.67	1.74	1.50	1.60	1.67
DuraVision® Plus Mirror UV	5450	5900	-	-	5600	6050	-
DuraVision® Plus Sun UV	5350	5800	6500	-	5500	5950	-
DuraVision® Plus Flash UV	5450	5900	6600	-	5600	6050	-
DuraVision® Plus Gold UV	5550	6000	6700	-	5700	6150	-
DuraVision® Plus Platinum UV	5350	5800	6500	-	5500	5950	-
Polarized	1.50	1.60	1.67	1.74			
DuraVision® Plus Mirror UV	5250	5700	-	-			
DuraVision® Plus Sun UV	5150	5600	6300	-			
DuraVision® Plus Flash UV	5250	5700	6400	-			
DuraVision® Plus Gold UV	5350	5800	6500	-			
DuraVision® Plus Platinum UV	5150	5600	6300	-			

* Require customer's date of birth (e.g. 14 Feb 1970), Frame Shape and Size, Pantoscopic angle, Vertex Distance, Bow Angle, FrameFit value (FFV from -1.0 to 6.0) for lens optimization.

* Require customer's Vision Profile: Balance, Intermediate or Near for lifestyle optimization.

* AdaptiveSun Blue is only available in 1.5 and 1.6 index.

* AdaptiveSun and AdaptiveSun Polarized is not available in DuraVision Mirror UV Strong Blue.

* Tinted colour is required for DuraVision Mirror UV, DuraVision Sun UV, DuraVision Flash UV (except for PhotoFusion X)

BlueGuard / UVProtect	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+6.00 to -7.00	-6.00	-7.00	+0.75 to +3.50	75	3
1.60	+10.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -17.00	-6.00	-17.00	+0.75 to +3.50	70	3
1.74	+16.00 to -20.00	-6.00	-20.00	+0.75 to +3.50	70	3
PhotoFusion X	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	75	3
1.74	+9.00 to -14.00	-6.00	-14.00	+0.75 to +3.50	75	3
AdaptiveSun	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	70	3
AdaptiveSun Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.75 to +3.50	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.75 to +3.50	80	3
Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.75 to +3.50	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.75 to +3.50	80	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	80	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Price quoted is as per pair in Ringgit Malaysia and subject to prevailing tax. Price is subject to change without prior notice.

ZEISS Single Vision Individual Sport	Now with BlueGuard					
	BlueGuard			PhotoFusion X		
	1.50	1.60	1.67	1.50	1.60	1.67
DuraVision® Plus Mirror UV	1500	1750	-	-	-	-
DuraVision® Plus Sun UV	1400	1650	1850	-	-	-
DuraVision® Plus Flash UV	1500	1750	1950	1850	2070	2310
DuraVision® Plus Gold UV	1500	1750	1950	1950	2170	2410
DuraVision® Plus Platinum UV	1300	1550	1750	1750	1970	2210

** List price for DuraVision Flash, DuraVision Sun, DuraVision Mirror on UV Protect is inclusive of Standard Tint colour, excluding Syklet and Special Tint.*

** Top-up RM100/pair for Skylet (Fun, Road and Sport) & Special Tint (ProGolf, ProGolf Gradient and Double Gradient Tints)*

*** PhotoFusion X 1.74 available in Grey, Extra Grey, Brown and Black only.*

AdaptiveSun	AdaptiveSun			AdaptiveSun Polarized		
	1.50	1.60	1.67	1.50	1.60	1.67
DuraVision® Plus Mirror UV	2050	2270	-	2200	2400	-
DuraVision® Plus Sun UV	1950	2170	2410	2100	2300	-
DuraVision® Plus Flash UV	2050	2270	2510	2200	2400	-
DuraVision® Plus Gold UV	2150	2370	2610	2300	2500	-
DuraVision® Plus Platinum UV	1950	2170	2410	2100	2300	-
Polarized	1.50	1.60	1.67			
DuraVision® Plus Mirror UV	1850	2070	-			
DuraVision® Plus Sun UV	1750	1970	2210			
DuraVision® Plus Flash UV	1850	2070	2310			
DuraVision® Plus Gold UV	1950	2170	2410			
DuraVision® Plus Platinum UV	1750	1970	2210			

** Require Frame Shape and Size, Pantoscopic angle, Bow Angle, lens base curve (optional) for lens optimization.*

** AdaptiveSun Blue is only available in 1.5 and 1.6 index.*

** AdaptiveSun and AdaptiveSun Polarized is not available in DuraVision Mirror UV Strong Blue.*

** Tinted colour is required for DuraVision Mirror UV, DuraVision Sun UV, DuraVision Flash UV (except for PhotoFusion X)*

UVProtect / PhotoFusion X / AdaptiveSun / AdaptiveSun Polarized / Polarized	SPH	CYL	COM RX.	DIA	PRISM
1.50	+5.00 to -5.00	-4.00	-5.00	75/80	3
1.60	+6.00 to -6.00	-6.00	-6.00	75/80	3
1.67	+6.00 to -6.00	-6.00	-6.00	75/80	3

** Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.*

** You are advised to use VISUSTORE® to verify the availability of maximum diameter.*

ZEISS Digital Lens Individual Sport	UVProtect			PhotoFusion X		
	1.50	1.60	1.67	1.50	1.60	1.67
DuraVision® Plus Mirror UV	1600	1850	-	-	-	-
DuraVision® Plus Sun UV	1500	1750	2150	-	-	-
DuraVision® Plus Flash UV	1600	1850	2250	1890	2250	2550
DuraVision® Plus Gold UV	1600	1850	2250	1990	2350	2650
DuraVision® Plus Platinum UV	1400	1650	2050	1790	2150	2450

* List price for DuraVision Flash, DuraVision Sun, DuraVision Mirror on UV Protect is inclusive of Standard Tint colour, excluding Syklet and Special Tint.
 * Top-up RM100/pair for Skylet (Fun, Road and Sport) & Special Tint (ProGolf, ProGolf Gradient and Double Gradient Tints)

AdaptiveSun	AdaptiveSun			AdaptiveSun Polarized		
	1.50	1.60	1.67	1.50	1.60	1.67
DuraVision® Plus Mirror UV	2090	2450	-	2240	2600	-
DuraVision® Plus Sun UV	1990	2350	2650	2140	2500	-
DuraVision® Plus Flash UV	2090	2450	2750	2240	2600	-
DuraVision® Plus Gold UV	2190	2550	2850	2340	2700	-
DuraVision® Plus Platinum UV	1990	2350	2650	2140	2500	-
Polarized	1.50	1.60	1.67			
DuraVision® Plus Mirror UV	1890	2250	-			
DuraVision® Plus Sun UV	1790	2150	2450			
DuraVision® Plus Flash UV	1890	2250	2550			
DuraVision® Plus Gold UV	1990	2350	2650			
DuraVision® Plus Platinum UV	1790	2150	2450			

* Require Frame Shape and Size, Pantoscopic angle, Bow Angle, Back Vertex Distance, lens base curve (optional) for lens optimization.
 * AdaptiveSun Blue is only available in 1.5 and 1.6 index.
 * AdaptiveSun and AdaptiveSun Polarized is not available in DuraVision Mirror UV Strong Blue.
 * Tinted colour is required for DuraVision Mirror UV, DuraVision Sun UV, DuraVision Flash UV (except for PhotoFusion X)

UVProtect / PhotoFusion X / AdaptiveSun / AdaptiveSun Polarized / Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -5.00	-4.00	-5.00	+0.50 to +1.25	75/80	3
1.60	+6.00 to -6.00	-6.00	-6.00	+0.50 to +1.25	75/80	3
1.67	+6.00 to -6.00	-6.00	-6.00	+0.50 to +1.25	75/80	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.
 * You are advised to use VISUSTORE® to verify the availability of maximum diameter.

ZEISS Progressive Individual Sport	Now with BlueGuard					
	UVProtect			PhotoFusion X		
	1.50	1.60	1.67	1.50	1.60	1.67
DuraVision® Plus Mirror UV	4590	5070	-	-	-	-
DuraVision® Plus Sun UV	4490	4970	5450	-	-	-
DuraVision® Plus Flash UV	4590	5070	5550	4950	5470	5900
DuraVision® Plus Gold UV	4590	5070	5550	5050	5570	6000
DuraVision® Plus Platinum UV	4390	4870	5350	4850	5370	5800

* List price for DuraVision Flash, DuraVision Sun, DuraVision Mirror on UV Protect is inclusive of Standard Tint colour, excluding Syklet and Special Tint.

* Top-up RM100/pair for Skylet (Fun, Road and Sport) & Special Tint (ProGolf, ProGolf Gradient and Double Gradient Tints)

AdaptiveSun	AdaptiveSun			AdaptiveSun Polarized		
	1.50	1.60	1.67	1.50	1.60	1.67
DuraVision® Plus Mirror UV	5150	5670	-	5300	5820	-
DuraVision® Plus Sun UV	5050	5570	6000	5200	5720	-
DuraVision® Plus Flash UV	5150	5670	6100	5300	5820	-
DuraVision® Plus Gold UV	5250	5770	6200	5400	5920	-
DuraVision® Plus Platinum UV	5050	5570	6000	5200	5720	-
Polarized	1.50	1.60	1.67			
DuraVision® Plus Mirror UV	4950	5470	-			
DuraVision® Plus Sun UV	4850	5370	5800			
DuraVision® Plus Flash UV	4950	5470	5900			
DuraVision® Plus Gold UV	5050	5570	6000			
DuraVision® Plus Platinum UV	4850	5370	5800			

* Require Frame Shape and Size, Pantoscopic angle, Bow Angle, FrameFit value (FFV:0-6), lens base curve (optional) for lens optimization.

* AdaptiveSun Blue is only available in 1.5 and 1.6 index.

* AdaptiveSun and AdaptiveSun Polarized is not available in DuraVision Mirror UV Strong Blue.

* Tinted colour is required for DuraVision Mirror UV, DuraVision Sun UV, DuraVision Flash UV (except for PhotoFusion X)

UVProtect / PhotoFusion X / AdaptiveSun / AdaptiveSun Polarized / Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -5.00	-4.00	-5.00	+0.75 to +3.50	75/80	3
1.60	+6.00 to -6.00	-6.00	-6.00	+0.75 to +3.50	75/80	3
1.67	+6.00 to -6.00	-6.00	-6.00	+0.75 to +3.50	75/80	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Now with BlueGuard								
ZEISS Single Vision DriveSafe Individual	UVProtect				PhotoFusion X			
	1.50	1.60	1.67	1.74	1.50	1.60	1.67	1.74**
DuraVision® Plus DriveSafe UV	1330	1670	2000	3000	2300	2980	3660	3900
AdaptiveSun	AdaptiveSun				AdaptiveSun Polarized			
	1.50	1.60	1.67	1.74	1.50	1.60	1.67	
DuraVision® Plus Sun UV	2420	3100	3710	-	2570	3250	-	
DuraVision® Plus Gold UV	2620	3300	3910	-	2770	3450	-	
DuraVision® Plus Platinum UV	2420	3100	3710	-	2570	3250	-	
Polarized	1.50	1.60	1.67	1.74				
DuraVision® Plus Sun UV	2300	2980	3660	-				
DuraVision® Plus Gold UV	2500	3180	3860	-				
DuraVision® Plus Platinum UV	2300	2980	3660	-				

* Require Frame Shape and Size, Pantoscopic angle, Vertex Distance, Bow Angle for lens optimization.

* AdaptiveSun Blue is not available for DriveSafe products.

**1.74 is available in PhotoFusion X (Black, Brown, Grey, Extra Grey) only

UVProtect	SPH	CYL	COM RX.	DIA	PRISM
1.50	+6.00 to -6.00	-4.00	-6.00	75	3
1.60	+10.00 to -10.00	-6.00	-10.00	75	3
1.67	+10.00 to -17.00	-6.00	-17.00	70	3
1.74	+16.00 to -20.00	-6.00	-20.00	70	3
PhotoFusion X	SPH	CYL	COM RX.	DIA	PRISM
1.50	+6.00 to -6.00	-4.00	-6.00	75	3
1.60	+8.00 to -10.00	-6.00	-10.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	75	3
1.74	+9.00 to -14.00	-6.00	-14.00	75	3
AdaptiveSun	SPH	CYL	COM RX.	DIA	PRISM
1.50	+6.00 to -6.00	-4.00	-6.00	75	3
1.60	+8.00 to -10.00	-6.00	-10.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	70	3
AdaptiveSun Polarized	SPH	CYL	COM RX.	DIA	PRISM
1.50	+8.00 to -8.00	-4.00	-8.00	80	3
1.60	+10.00 to -11.00	-6.00	-11.00	80	3
Polarized	SPH	CYL	COM RX.	DIA	PRISM
1.50	+8.00 to -8.00	-4.00	-8.00	80	3
1.60	+10.00 to -11.00	-6.00	-11.00	80	3
1.67	+11.00 to -12.00	-6.00	-12.00	80	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Now with BlueGuard						
ZEISS Single Vision DriveSafe Individual Sport	UVProtect			PhotoFusion X		
	1.50	1.60	1.67	1.50	1.60	1.67
DuraVision® Plus DriveSafe UV	1260	1580	1900	2180	2820	3460
AdaptiveSun	AdaptiveSun			AdaptiveSun Polarized		
	1.50	1.60	1.67	1.50	1.60	1.67
DuraVision® Plus Sun UV	2290	2980	3590	2440	3130	-
DuraVision® Plus Gold UV	2490	3180	3790	2640	3330	-
DuraVision® Plus Platinum UV	2290	2980	3590	2440	3130	-
Polarized	1.50	1.60	1.67			
DuraVision® Plus Sun UV	2180	2820	3460			
DuraVision® Plus Gold UV	2380	3020	3660			
DuraVision® Plus Platinum UV	2180	2820	3460			

* Require Frame Shape and Size, Pantoscopic angle, Bow Angle, lens base curve (optional) for lens optimization.

* AdaptiveSun Blue is not available for DriveSafe product.

* Tinted colour on UVProtect range is not available.

UVProtect / PhotoFusion X / AdaptiveSun / AdaptiveSun Polarized / Polarized	SPH	CYL	COM RX.	DIA	PRISM
1.50	+5.00 to -5.00	-4.00	-5.00	75/80	3
1.60	+6.00 to -6.00	-6.00	-6.00	75/80	3
1.67	+6.00 to -6.00	-6.00	-6.00	75/80	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Now with BlueGuard								
ZEISS Progressive DriveSafe Individual	UVProtect				PhotoFusion X			
	1.50	1.60	1.67	1.74	1.50	1.60	1.67	1.74**
DuraVision® Plus DriveSafe UV	3030	3520	4060	4570	3530	4020	4560	5070
AdaptiveSun	AdaptiveSun				AdaptiveSun Polarized			
	1.50	1.60	1.67	1.74	1.50	1.60	1.67	
DuraVision® Plus Sun UV	3600	4320	4990	-	3750	4470	-	
DuraVision® Plus Gold UV	3800	4520	5190	-	3950	4670	-	
DuraVision® Plus Platinum UV	3600	4320	4990	-	3750	4470	-	
Polarized	1.50	1.60	1.67	1.74				
DuraVision® Plus Sun UV	3400	4120	4790	-				
DuraVision® Plus Gold UV	3600	4320	4990	-				
DuraVision® Plus Platinum UV	3400	4120	4790	-				

* Require Frame Shape and Size, Pantoscopic angle, Vertex Distance, Bow Angle, FrameFit Value (FFV: 0.0 - 6.0) for lens optimization.

* AdaptiveSun Blue is not available for DriveSafe products.

**1.74 is available in PhotoFusion X (Black, Brown, Grey, Extra Grey) only

UVProtect	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+6.00 to -7.00	-6.00	-7.00	+0.75 to +3.50	75	3
1.60	+10.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -17.00	-6.00	-17.00	+0.75 to +3.50	70	3
1.74	+16.00 to -20.00	-6.00	-20.00	+0.75 to +3.50	70	3
PhotoFusion X	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	75	3
1.74	+9.00 to -14.00	-6.00	-14.00	+0.75 to +3.50	75	3
AdaptiveSun	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	70	3
AdaptiveSun Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.75 to +3.50	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.75 to +3.50	80	3
Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.75 to +3.50	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.75 to +3.50	80	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	80	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Now with BlueGuard						
ZEISS Progressive DriveSafe Individual Sport	UVProtect			PhotoFusion X		
	1.50	1.60	1.67	1.50	1.60	1.67
DuraVision® Plus DriveSafe UV	2990	3460	3940	3310	3980	4620
AdaptiveSun	AdaptiveSun			AdaptiveSun Polarized		
	1.50	1.60	1.67	1.50	1.60	1.67
DuraVision® Plus Sun UV	3410	4080	4720	3560	4230	-
DuraVision® Plus Gold UV	3610	4280	4920	3760	4430	-
DuraVision® Plus Platinum UV	3410	4080	4720	3560	4230	-
Polarized	1.50	1.60	1.67			
DuraVision® Plus Sun UV	3310	3980	4620			
DuraVision® Plus Gold UV	3510	4180	4820			
DuraVision® Plus Platinum UV	3310	3980	4620			

* Require Frame Shape and Size, Pantoscopic angle, Bow Angle, FrameFit Value (FFV:0–6), lens base curve (optional) for lens optimization.

* AdaptiveSun Blue is not available for DriveSafe product.

* Tinted colour on UVProtect range is not available.

UVProtect / PhotoFusion X / AdaptiveSun / AdaptiveSun Polarized / Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -5.00	-4.00	-5.00	+0.75 to +3.50	75/80	3
1.60	+6.00 to -6.00	-6.00	-6.00	+0.75 to +3.50	75/80	3
1.67	+6.00 to -6.00	-6.00	-6.00	+0.75 to +3.50	75/80	3


* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

* You are advised to use VISUSTORE® to verify the availability of maximum diameter.



// OFFICE LENS INDIVIDUAL

Maximum Intermediate Distance technology by ZEISS - at a glance.



ZEISS OfficeLenses Individual	BlueGuard			
	1.50	1.60	1.67	1.74
DuraVision® Plus Gold UV	1620	1950	2320	2950
DuraVision® Plus Platinum UV	1420	1750	2120	2750
	UVProtect			
	1.50	1.60	1.67	1.74
DuraVision® Plus Gold UV	1520	1850	2220	2850
DuraVision® Plus Platinum UV	1320	1650	2020	2650

* Require Frame Shape and Size, FrameFit value (FFV from 0.0 to 6.0) and Maximum Intermediate Distance (MID) for lens optimization
* MID: Book (100cm), Near (200cm) or Room (400cm). Recommended MID is 100 or 200.

BlueGuard	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	70	3
1.74	+9.00 to -14.00	-6.00	-14.00	+0.75 to +3.50	70	3
UVProtect	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	70	3
1.74	+9.00 to -14.00	-6.00	-14.00	+0.75 to +3.50	70	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.
* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Price quoted is as per pair in Ringgit Malaysia and subject to prevailing tax. Price is subject to change without prior notice.



// INDIVIDUAL MYOPIA MANAGEMENT

ZEISS Myovision® Pro Individual	Now with BlueGuard					
	UVProtect			PhotoFusion X		
	1.50	1.60	1.67	1.50	1.60	1.67
DuraVision® Plus BlueProtect UV	1220	1460	1820	-	-	-
DuraVision® Plus Platinum UV	1120	1360	1720	1600	1840	2200

* Require Frame Shape and Size, Pantoscopic angle, Vertex Distance, Bow Angle for lens optimization.
* Tinted colour is unavailable.

UVProtect	SPH	CYL	COM RX.	DIA	PRISM
1.50	0.00 to -6.00	-2.00	-6.00	70	3
1.60	0.00 to -9.50	-4.00	-9.50	70	3
1.67	0.00 to -9.50	-4.00	-9.50	70	3
PhotoFusion X	SPH	CYL	COM RX.	DIA	PRISM
1.50	0.00 to -6.00	-2.00	-6.00	70	3
1.60	0.00 to -9.50	-4.00	-9.50	70	3
1.67	0.00 to -9.50	-4.00	-9.50	70	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.
* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

ZEISS Myokids® Pro Individual	Now with BlueGuard					
	UVProtect			PhotoFusion X		
	1.50	1.60	1.67	1.50	1.60	1.67
DuraVision® Plus BlueProtect UV	1280	1520	1880	-	-	-
DuraVision® Plus Platinum UV	1180	1420	1780	1660	1900	2260

* Require Frame Shape and Size, Pantoscopic angle, Vertex Distance, Bow Angle for lens optimization.
* PhotoFusion Extra Grey, Pioneer and PhotoFusion Blue require 14 days delivery.
* Tinted colour is unavailable.

UVProtect	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	0.00 to -7.00	-4.00	-7.00	+1.00 to +2.50	75	3
1.60	0.00 to -10.00	-6.00	-10.00	+1.00 to +2.50	75	3
1.67	0.00 to -12.00	-6.00	-12.00	+1.00 to +2.50	70	3
PhotoFusion X	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	0.00 to -7.00	-4.00	-7.00	+1.00 to +2.50	75	3
1.60	0.00 to -10.00	-6.00	-10.00	+1.00 to +2.50	75	3
1.67	0.00 to -12.00	-6.00	-12.00	+1.00 to +2.50	70	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.
* You are advised to use VISUSTORE® to verify the availability of maximum diameter.



// INDIVIDUAL MINERAL LENS

ZEISS SmartLife Progressive Individual 3 Mineral Lens	Clear		Umbramatic
	1.60	1.80	1.60
Super ET	5100	6400	5700

* Require customer's date of birth (e.g. 14 Feb 1970), Frame Shape and Size, Pantoscopic angle, Vertex Distance, Bow Angle, FrameFit value (FFV from -1.0 to 6.0) for lens optimization.
* Require customer's Vision Profile: Balance, Intermediate or Near for lifestyle optimization.
* Require 14 days delivery.
* Tinted colour is unavailable.

Clear	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.60	+6.50 to -11.00	-6.00	-11.00	+0.75 to +3.50	70	3
1.80	0.00 to -20.00	-8.00	-20.00	+0.75 to +3.50	65	3
Umbramatic	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.60	+6.50 to -9.00	-4.00	-9.00	+0.75 to +3.50	70	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.
You are advised to use VISUSTORE® to verify the availability of maximum diameter.




ZEISS SmartLife Digital Individual 3 Mineral Lens	Clear		Umbramatic
	1.60	1.80	1.60
Super ET	2000	2800	2600

* Require customer's date of birth (e.g. 14 Feb 1970), Frame Shape and Size, Pantoscopic angle, Vertex Distance, Bow Angle, FrameFit value (FFV from -1.0 to 6.0) for lens optimization.
* Require customer's date of birth (e.g. 14 Feb 1970), Frame Shape and Size for lens optimization.
* Require 14 days delivery.
* Tinted colour is unavailable.

Clear	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.60	+6.50 to -11.00	-6.00	-11.00	+0.50 to +1.25	70	3
1.80	0.00 to -20.00	-8.00	-20.00	+0.50 to +1.25	65	3
Umbramatic	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.60	+6.50 to -9.00	-4.00	-9.00	+0.50 to +1.25	70	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.
You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Stamp and Engravings



ZEISS SmartLife Lenses

Lens Engraving

Nasal Engraving	SmartLife	SmartLife PRO
ZEISS Single Vision SmartLife	SL	VSP
ZEISS Single Vision SmartLife Asiana	SLA	
ZEISS Single Vision SmartLife Individual 3	VI3	VP3
ZEISS SmartLife Digital	DSL	DSP
ZEISS SmartLife Digital Asiana	DLA	
ZEISS Digital SmartLife Individual 3	DI3	DP3
ZEISS Progressive SmartLife Pure	SPU	SPR
ZEISS Progressive SmartLife Plus	SPL	SPP
ZEISS Progressive SmartLife Plus Asiana	SPA	
ZEISS Progressive SmartLife Superb	SSB	SPS
ZEISS Progressive SmartLife Individual 3	SI3	SP3