

Seeing beyond

ZEISS

ZEISS ClearMind Individual 3

Experience truly personalised vision correction with the Individual 3 customisation level — crafted exclusively for your unique perspective.

See the world your way.

Lenses customized especially for your eyes.

Unlock the full potential of your vision with the **Individual 3 customization level** for ZEISS ClearMind lenses. We use ZEISS Intelligence Augmented Design technology to enhance your personal, anatomical and physiological data with **12.5 million data points** about real-life, age-related visual behaviour collected from by ZEISS Global Vision Study.*

Thoughtfully designed to suit your unique anatomy, specific vision needs and age-related visual behaviour, these lenses adapt perfectly to your favourite frames to deliver extremely clear vision and comfort.

Our partner eye care professionals use the latest ZEISS measuring technology to personalise the perfect ZEISS ClearMind lenses for you and your eyes.



* ZEISS Global Vision Study 2020/21.

** for intermediate distance (with ZEISS Progressive ClearMind Individual 3 compared to ZEISS Progressive ClearMind Superb Lenses). Calculation & Analysis by Technology & Innovation, Carl Zeiss Vision International GmbH, DE 2025 (unpublished, data on file)

A deep look in *your* eyes.

Ready to make your ZEISS ClearMind lenses even better?

Dedicated to providing you with the best possible visual experience, we use advanced ZEISS Individual Luminance Design technology to create a personalized visual optical profile. This takes the unique characteristics of your eyes into account, including the size of your pupils.

We use this data to optimize your prescription and upgrade your lenses to






ZEISS ClearMind PRO lenses.

Ask your eye care expert about this premium option.



ZEISS Digital & Progressive ClearMind lenses

Tailored solutions to serve specific needs

DIGITAL	35+ For individuals 35+ who experience blurry vision, near or far, along with increased eye fatigue due to age-related changes. Accommodative support ranging from +0.50D to 1.25D.	Lens options	Vision zones	Frame choice
		OUR BEST SOLUTION Individual 3 tier 	• Larger clear zones when fitted to individual anatomy • Accommodative support in the lower lens part for relaxed reading Smooth peripheral blur zones • Optimised for digital device use • Optimised for individuals wearing addition power lenses for the first time • Optimised for support of the cognitive load reduction	 For any frame of your choice, adapted to your individual anatomy and special vision needs.
		OUR BETTER SOLUTION Superb tier 	• Large clear zones • Accommodative support in the lower lens part for relaxed reading Smooth peripheral blur zones • Optimised for digital device use • Optimized for individuals wearing addition power lenses for the first time • Optimised for support of the cognitive load reduction	 For any frame of your choice, adapted to your individual anatomy.
PROGRESSIVE	45+ For individuals 45+ who need seamless vision from near to far and at all distances. From +0.75D up to +4.00D addition power.	OUR BEST SOLUTION Individual 3 tier 	• Larger clear zones than ZEISS • Progressive ClearMind Superb when fitted to individual anatomy & vision needs. Smooth peripheral blur zones. • Optimised for digital device use. • Optimised for support of the cognitive load reduction. Up to 19% larger vision zones.*	 For any frame of your choice, adapted to your individual anatomy and special vision needs.
		OUR BETTER SOLUTION Superb tier 	• Larger clear zones than ZEISS • Progressive ClearMind Plus when fitted to individual anatomy. • Smooth peripheral blur. • Optimised for digital device use. • Optimised for support of the cognitive load reduction.	 For any frame of your choice, adapted to your individual anatomy.
		OUR GOOD SOLUTION Plus tier 	• Larger clear zones than ZEISS • Progressive ClearMind Pure when fitted in near vision zone. • Smooth peripheral blur. • Optimised for digital device use. • Optimised for support of the cognitive load reduction.	 For any frame with variable fitting heights up from 13 mm.
		OUR BASIC SOLUTION Pure tier 	• Large clear zones. • Smooth peripheral blur. • Optimised for digital device usage. • Optimised for support of the cognitive load reduction.	 For standard frames available in 3 fitting heights (14/16/18 mm).

* For intermediate distance (with ZEISS Progressive ClearMind Individual 3 compared to ZEISS Progressive ClearMind Superb Lenses). Calculation & Analysis by Technology & Innovation, Carl Zeiss Vision International GmbH, DE 2025 (unpublished, data on file)

ZEISS ClearMind Single Vision Individual 3	BlueGuard							
	1.50	1.60	1.67	1.74				
DuraVision® Plus Gold UV	1980	2230	2470	3250				
DuraVision® Plus Platinum UV	1780	2030	2270	3050	Now with BlueGuard			
Clear PhotoFusion X Tinted	UVProtect				PhotoFusion X			
	1.50	1.60	1.67	1.74	1.50	1.60	1.67	1.74**
DuraVision® Plus Mirror UV	1880	2130	-	-	-	-	-	-
DuraVision® Plus Sun UV	1780	2030	2270	-	-	-	-	-
DuraVision® Plus Flash UV	1880	2130	2370	-	2380	2630	2870	-
DuraVision® Plus Gold UV	1880	2130	2370	3150	2480	2730	2970	3750
DuraVision® Plus Platinum UV	1680	1930	2170	2950	2280	2530	2770	3550

* List price for DuraVision® Plus Flash, DuraVision® Plus Sun, DuraVision® Plus Mirror on UV Protect is inclusive of Standard Tint colour, excluding Syklet and Special Tint.

* Top-up RM100/pair for Syklet (Fun, Road and Sport) & Special Tint (ProGolf, ProGolf Gradient and Double Gradient Tints)

**PhotoFusion X 1.74 available in Grey, Extra Grey, Brown and Black only.

AdaptiveSun	AdaptiveSun				AdaptiveSun Polarized		
	1.50	1.60	1.67	1.74	1.50	1.60	1.67
DuraVision® Plus Mirror UV	2580	2830	-	-	2730	2980	-
DuraVision® Plus Sun UV	2480	2730	2970	-	2630	2880	-
DuraVision® Plus Flash UV	2580	2830	3070	-	2730	2980	-
DuraVision® Plus Gold UV	2680	2930	3170	-	2830	3080	-
DuraVision® Plus Platinum UV	2480	2730	2970	-	2630	2880	-
Polarized	1.50	1.60	1.67	1.74			
DuraVision® Plus Mirror UV	2380	2630	-	-			
DuraVision® Plus Sun UV	2280	2530	2770	-			
DuraVision® Plus Flash UV	2380	2630	2870	-			
DuraVision® Plus Gold UV	2480	2730	2970	-			
DuraVision® Plus Platinum UV	2280	2530	2770	-			

* Require customer's date of birth (e.g. 14 Feb 1996). Frame Shape and Size, Pantoscopic angle, Wrap angle and Back Vertex Distance, for lens optimization

* AdaptiveSun Blue is only available in 1.5 and 1.6 index.

* AdaptiveSun and AdaptiveSun Polarized is not available in DuraVision® Plus Mirror UV Strong Blue.

* Tinted colour is required for DuraVision® Plus Mirror UV, DuraVision® Plus Sun UV, DuraVision® Plus Flash UV (except for PhotoFusion X)

BlueGuard / UVProtect	SPH	CYL	COM RX.	DIA	PRISM
1.50	+6.00 to -6.00	-4.00	-6.00	75	3
1.60	+10.00 to -10.00	-6.00	-10.00	75	3
1.67	+10.00 to -17.00	-6.00	-17.00	70	3
1.74	+16.00 to -20.00	-6.00	-20.00	70	3
PhotoFusion X	SPH	CYL	COM RX.	DIA	PRISM
1.50	+6.00 to -6.00	-4.00	-6.00	75	3
1.60	+8.00 to -10.00	-6.00	-10.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	70	3
1.74	+9.00 to -14.00	-6.00	-14.00	75	3
AdaptiveSun	SPH	CYL	COM RX.	DIA	PRISM
1.50	+6.00 to -6.00	-4.00	-6.00	75	3
1.60	+8.00 to -10.00	-6.00	-10.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	70	3
AdaptiveSun Polarized	SPH	CYL	COM RX.	DIA	PRISM
1.50	+6.00 to -6.00	-4.00	-6.00	80	3
1.60	+10.00 to -11.00	-6.00	-11.00	80	3
Polarized	SPH	CYL	COM RX.	DIA	PRISM
1.50	+8.00 to -8.00	-4.00	-8.00	80	3
1.60	+10.00 to -11.00	-6.00	-11.00	80	3
1.67	+11.00 to -12.00	-6.00	-12.00	80	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Price quoted is as per pair in Ringgit Malaysia and subject to prevailing tax. Price is subject to change without prior notice.

ZEISS ClearMind Digital Individual 3	BlueGuard							
	1.50	1.60	1.67	1.74				
DuraVision® Plus Gold UV	2150	2390	2800	3270				
DuraVision® Plus Platinum UV	1950	2190	2600	3070	Now with BlueGuard			
Clear PhotoFusion X Tinted	UVProtect				PhotoFusion X			
	1.50	1.60	1.67	1.74	1.50	1.60	1.67	1.74**
DuraVision® Plus Mirror UV	2050	2290	-	-	-	-	-	-
DuraVision® Plus Sun UV	1950	2190	2600	-	-	-	-	-
DuraVision® Plus Flash UV	2050	2290	2700	-	2550	2790	3200	-
DuraVision® Plus Gold UV	2050	2290	2700	3170	2650	2890	3300	3770
DuraVision® Plus Platinum UV	1850	2090	2500	2970	2450	2690	3100	3570

* List price for DuraVision® Plus Flash, DuraVision® Plus Sun, DuraVision® Plus Mirror on UV Protect is inclusive of Standard Tint colour, excluding Syklet and Special Tint.

* Top-up RM100/pair for Syklet (Fun, Road and Sport) & Special Tint (ProGolf, ProGolf Gradient and Double Gradient Tints)

** PhotoFusion X 1.74 available in Grey, Extra Grey, Brown and Black only.

AdaptiveSun	AdaptiveSun				AdaptiveSun Polarized		
	1.50	1.60	1.67	1.74	1.50	1.60	1.67
DuraVision® Plus Mirror UV	2750	2990	-	-	2900	3140	-
DuraVision® Plus Sun UV	2650	2890	3300	-	2800	3040	-
DuraVision® Plus Flash UV	2750	2990	3400	-	2900	3140	-
DuraVision® Plus Gold UV	2850	3090	3500	-	3000	3240	-
DuraVision® Plus Platinum UV	2650	2890	3300	-	2800	3040	-
Polarized	1.50	1.60	1.67	1.74			
DuraVision® Plus Mirror UV	2550	2790	-	-			
DuraVision® Plus Sun UV	2450	2690	3100	-			
DuraVision® Plus Flash UV	2550	2790	3200	-			
DuraVision® Plus Gold UV	2650	2890	3300	-			
DuraVision® Plus Platinum UV	2450	2690	3100	-			

* Require customer's date of birth (eg. 14 Feb 1996), Frame Shape and Size, Pantoscopic angle for lens optimization, FFV

* AdaptiveSun Blue is only available in 1.5 and 1.6 index.

* AdaptiveSun and AdaptiveSun Polarized is not available in DuraVision® Plus Mirror UV Strong Blue.

* Tinted colour is required for DuraVision® Plus Mirror UV, DuraVision® Plus Sun UV, DuraVision® Plus Flash UV (except for PhotoFusion X)

BlueGuard / UVProtect	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+6.00 to -7.00	-6.00	-7.00	+0.50 to +1.25	75	3
1.60	+10.00 to -10.00	-6.00	-10.00	+0.50 to +1.25	75	3
1.67	+8.00 to -17.00	-6.00	-17.00	+0.50 to +1.25	70	3
1.74	+16.00 to -20.00	-6.00	-20.00	+0.50 to +1.25	70	3
PhotoFusion X	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.50 to +1.25	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.50 to +1.25	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.50 to +1.25	75	3
1.74	+9.00 to -14.00	-6.00	-14.00	+0.50 to +1.25	75	3
AdaptiveSun	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.50 to +1.25	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.50 to +1.25	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.50 to +1.25	70	3
AdaptiveSun Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.50 to +1.25	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.50 to +1.25	80	3
Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.50 to +1.25	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.50 to +1.25	80	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.50 to +1.25	80	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Price quoted is as per pair in Ringgit Malaysia and subject to prevailing tax. Price is subject to change without prior notice.

ZEISS ClearMind Progressive Individual 3	BlueGuard							
	1.50	1.60	1.67	1.74				
DuraVision® Plus Gold UV	5238	5788	6538	7238				
DuraVision® Plus Platinum UV	5050	5600	6350	7050	Now with BlueGuard			
Clear PhotoFusion X Tinted	UVProtect				PhotoFusion X			
	1.50	1.60	1.67	1.74	1.50	1.60	1.67	1.74**
DuraVision® Plus Mirror UV	5150	5700	-	-	-	-	-	-
DuraVision® Plus Sun UV	5050	5600	6350	-	-	-	-	-
DuraVision® Plus Flash UV	5150	5700	6450	-	5650	6200	6950	-
DuraVision® Plus Gold UV	5150	5700	6450	7150	5750	6300	7050	7750
DuraVision® Plus Platinum UV	4950	5500	6250	6950	5550	6100	6850	7550

* List price for DuraVision® Plus Flash, DuraVision® Plus Sun, DuraVision® Plus Mirror on UV Protect is inclusive of Standard Tint colour, excluding Syklet and Special Tint.

* Top-up RM100/pair for Skylet (Fun, Road and Sport) & Special Tint (ProGolf, ProGolf Gradient and Double Gradient Tints)

**PhotoFusion X 1.74 available in Grey, Extra Grey, Brown and Black only.

AdaptiveSun	AdaptiveSun				AdaptiveSun Polarized		
	1.50	1.60	1.67	1.74	1.50	1.60	1.67
DuraVision® Plus Mirror UV	5850	6400	-	-	6000	6550	-
DuraVision® Plus Sun UV	5750	6300	7050	-	5900	6450	-
DuraVision® Plus Flash UV	5850	6400	7150	-	6000	6550	-
DuraVision® Plus Gold UV	5950	6500	7250	-	6100	6650	-
DuraVision® Plus Platinum UV	5750	6300	7050	-	5900	6450	-
Polarized	1.50	1.60	1.67	1.74			
DuraVision® Plus Mirror UV	5650	6200	-	-			
DuraVision® Plus Sun UV	5550	6100	6850	-			
DuraVision® Plus Flash UV	5650	6200	6950	-			
DuraVision® Plus Gold UV	5750	6300	7050	-			
DuraVision® Plus Platinum UV	5550	6100	6850	-			

* Require customer's date of birth (e.g. 14 Feb 1970), Frame Shape and Size, Pantoscopic angle, Vertex Distance, Bow Angle, FrameFit value (FFV from -1.0 to 6.0) for lens optimization.

* Require customer's Vision Profile: Balance, Intermediate or Near for lifestyle optimization.

* AdaptiveSun Blue is only available in 1.5 and 1.6 index.

* AdaptiveSun and AdaptiveSun Polarized is not available in DuraVision® Plus Mirror UV Strong Blue.

* Tinted colour is required for DuraVision® Plus Mirror UV, DuraVision® Plus Sun UV, DuraVision® Plus Flash UV (except for PhotoFusion X)

BlueGuard / UVProtect	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+6.00 to -7.00	-6.00	-7.00	+0.75 to +3.50	75	3
1.60	+10.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -17.00	-6.00	-17.00	+0.75 to +3.50	70	3
1.74	+16.00 to -20.00	-6.00	-20.00	+0.75 to +3.50	70	3
PhotoFusion X	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	75	3
1.74	+9.00 to -14.00	-6.00	-14.00	+0.75 to +3.50	75	3
AdaptiveSun	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -7.00	-4.00	-7.00	+0.75 to +3.50	75	3
1.60	+6.00 to -10.00	-6.00	-10.00	+0.75 to +4.00	75	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	70	3
AdaptiveSun Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.75 to +3.50	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.75 to +3.50	80	3
Polarized	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.50	+5.00 to -8.00	-4.00	-8.00	+0.75 to +3.50	80	3
1.60	+6.00 to -11.00	-6.00	-11.00	+0.75 to +3.50	80	3
1.67	+8.00 to -12.00	-6.00	-12.00	+0.75 to +3.50	80	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

* You are advised to use VISUSTORE® to verify the availability of maximum diameter.

Price quoted is as per pair in Ringgit Malaysia and subject to prevailing tax. Price is subject to change without prior notice.

ZEISS ClearMind Progressive Individual 3 Mineral Lens	Clear		Umbramatic
	1.60	1.80	1.60
Super ET	5600	7050	6200

* Require customer's date of birth (e.g. 14 Feb 1970), Frame Shape and Size, Pantoscopic Angle, Back Vertex Distance, FrameFit (FFV: 0.0 - 6.0) for lens optimization.

* Require 14 days delivery.

* Tinted colour is unavailable.

Clear	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.60	+6.50 to -11.00	-6.00	-11.00	+0.75 to +3.50	70	3
1.80	0.00 to -20.00	-8.00	-20.00	+0.75 to +3.50	65	3
Umbramatic	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.60	+6.50 to -9.00	-4.00	-9.00	+0.75 to +3.50	70	3

* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

You are advised to use VISUSTORE® to verify the availability of maximum diameter.

ZEISS ClearMind Digital Individual 3 Mineral Lens	Clear		Umbramatic
	1.60	1.80	1.60
Super ET	2190	3070	2790

* Require customer's date of birth (e.g. 14 Feb 1970), Frame Shape and Size, Pantoscopic Angle, Back Vertex Distance, FrameFit (FFV: 0.0 - 6.0) for lens optimization.

* Require 14 days delivery.


* Tinted colour is unavailable.

Clear	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.60	+6.50 to -11.00	-6.00	-11.00	+0.50 to +1.25	70	3
1.80	0.00 to -20.00	-8.00	-20.00	+0.50 to +1.25	65	3
Umbramatic	SPH	CYL	COM RX.	ADD	DIA	PRISM
1.60	+6.50 to -9.00	-4.00	-9.00	+0.50 to +1.25	70	3


* Maximum diameter is subjected to combine RX. Maximum diameter typically is smaller with higher combined RX.

You are advised to use VISUSTORE® to verify the availability of maximum diameter.


Stamp and Engravings



ZEISS ClearMind Single Vision lenses




ZEISS Digital ClearMind Lenses lenses




ZEISS Digital ClearMind Individual 3 Lenses

Stamp and Engravings



ZEISS Progressive ClearMind Individual 3 lens



ZEISS ClearMind Lenses Lens Engraving

Nasal Engraving	ClearMind
ZEISS Single Vision ClearMind	VCMxx
ZEISS Single Vision ClearMind Individual 3	VCIdx
ZEISS Digital ClearMind	BCMxx
ZEISS Digital Asiana* ClearMind	DCAxx
ZEISS Digital ClearMind Individual 3	BCIdx
ZEISS Progressive ClearMind Pure	CPUxx
ZEISS Progressive ClearMind Plus	CMPxx
ZEISS Progressive ClearMind Plus Asiana*	CPAxx
ZEISS Progressive ClearMind Superb	CMSxx
ZEISS Progressive ClearMind Individual 3	CMIxx

*Asiana version only available for some countries.